Hepatitis C Testing	
Test	Qualitative assay for the detection of antibody to hepatitis C virus (anti-HCV) in
Description	human serum or plasma.
Test Use	As a screening assay to aid in the diagnosis of recent and/or past infection with hepatitis C virus (HCV).
Test	Virology
Department	Phone (860) 920-6662, FAX: (860) 920-6661
Methodology	Enzyme-linked Immunosorbent Assay (ELISA)
Availability	Test is performed 2 or 3 times weekly
Specimen	1 mL serum (preferred) or plasma derived from EDTA, lithium heparin, CPD,
Requirements	CP2D, CPDA-1, ACD or 4% citrate anticoagulants
Collection	To obtain collection kit, refer to Collection Kit Ordering Information.
kit/container	
Collection	Standard venipuncture technique
Instructions	
Specimen	Store specimens at room temperature (15-25 °C) or at 2-8° C for up to 7 days,
Handling &	including transit time. For longer storage of specimens, keep at -20° C for up to
Transport	4 weeks. Transport to laboratory with ice packs.
	Unlabeled specimen
Unacceptable	Specimens that have leaked or containers that have broken in transit
Conditions	Hemolyzed or heat-treated specimens
Requisition	Clinical Test Requisition OL-9B (select Hepatitis C Testing)
Form	
	Name and address of submitter. Two patient identifiers (ie.name, DOB, Acc.#,
Required	MRN), Town of residence (city, state, zip), specimen source/type, date
Information	collected, test(s) requested.
	Please ensure information on the requisition matches that on the specimen.
	The presence of anti-HCV does not constitute a diagnosis of hepatitis C disease
Limitations	and may be indicative of recent and/or past HCV infection.
	A nonreactive test result does not exclude the possibility of exposure to HCV.
	Levels of anti-HCV may be undetectable in early infection.
Additional	Nucleic acid amplification testing for Hepatitis C RNA is recommended for
Comments	patients with repeatedly reactive HCV antibody test results. Repeatedly reactive
	HCV antibody specimens are reflexed to Hepatitis C RNA testing when specimen
	volume is sufficient and specimen stability requirements are meet.

Revision: 1/2/2024