Fentanyl Plus ⁺ Drug Screen Revised 10-24-2023	
Test Description	Quantitative screening assay for the presumptive identification of opiates, opioids, prescription drugs, illicit drugs, and their metabolites in urine, including novel drug analogues and associated Tentatively Identified Compounds (TICs).
Test Use	To aid in the surveillance of illicit drug use in hospital emergency department admissions due to overdose.
Test Department	CT DPH Laboratory - Clinical Chemistry Phone: (860) 920-6592 Fax: (860) 920-6718
Methodology	Liquid Chromatography/Tandem Mass Spectrometry: LC MS/MS QTOF X500R
Availability	Daily, Monday-Friday
Specimen Requirements	Urine, 2.0-10.0 mL Low specimen volume may compromise repeat analysis or confirmatory testing.
Collection Kit/Container	Single-use plastic urine collection container, leak resistant. Containers and urine tube aliquots that contain additives/stabilizers cannot be tested.
Collection Instructions	Collect urine midstream (clean catch), by catheterization, or suprapubic aspiration. A sterile specimen is not required.
Specimen Handling & Transport	Store specimens at 2-8°C for up to 24 hours. If a delay is expected, store specimens at -20°C or lower. Transport specimens using ice packs.
Unacceptable Conditions	Unlabeled specimen. Insufficient specimen volume. Specimens that have leaked or containers that have broken in transit. Specimens not collected or handled as described above.
Requisition Form	Clinical Test Requisition OL-9B, under Test Not Listed (Specify), write in "Fentanyl Plus ⁺ Drug Screen".
Required Information	Name and address of submitter. Two patient identifiers (ie. name, DOB, Acc.#, MRN), name of provider, specimen type or source, date of collection, test requested. Please ensure patient information on the requisition matches that on the specimen.
Limitations	Screening is limited only to clinical specimens from patients that present at hospital emergency departments and meet the defined criteria for overdose.
Additional Comments	Screening is offered for surveillance purposes only under the CT DPH Injury Prevention Program and is not valid for medicolegal use. Screening alone may not identify novel compounds, with independent confirmation indicated.