

## Chlamydia & Gonorrhea Nucleic Acid Amplification Test

Revised 7/16/18

<b>Test description</b>	Qualitative detection and differentiation of ribosomal RNA (rRNA) from <i>Chlamydia trachomatis</i> and/or <i>Neisseria gonorrhoeae</i> in genital, extra genital and urine specimens.
<b>Test use</b>	As an aid to the diagnosis of chlamydial and gonococcal disease in symptomatic or asymptomatic individuals.
<b>Test Department</b>	Sexually Transmitted Diseases Molecular Diagnostics Phone: (860) 920-6696; Fax: (860) 920-6721
<b>Methodology</b>	Target amplification nucleic acid probe test (Aptima Combo 2 <sup>®</sup> Assay)
<b>Availability</b>	Daily, Monday - Friday
<b>Specimen Requirements</b>	Specimens must be received by the laboratory: <ul style="list-style-type: none"> <li>• Within 30 days of collection for urine, oropharyngeal swabs and rectal swabs.</li> <li>• Within 60 days of collection for vaginal, female endocervical and male urethral swabs.</li> </ul>
<b>Collection Kit/Container</b>	<ul style="list-style-type: none"> <li>• <u>Urine</u>: Aptima<sup>®</sup> Urine Specimen Collection Kit. <u>Vaginal, oropharyngeal and rectal swabs</u>: Aptima<sup>®</sup> Multitest Swab Specimen Collection Kit. <u>Female endocervical and male urethral swabs</u>: Aptima<sup>®</sup> Unisex Swab Specimen Collection Kit</li> </ul> <p>To obtain collection kits, refer to Collection Kit Ordering Information. Store collection kits at room temperature prior to use.</p>
<b>Collection instructions</b>	Follow instructions included with each specific collection kit. Specimen must be collected before the expiration date on the collection kit.
<b>Specimen Handling &amp; Transport</b>	Once collected, transport and store the specimens at 2°C-30°C (36°F -86°F). Avoid temperature extremes.
<b>Unacceptable Conditions</b>	Unlabeled specimens. Specimens that have leaked or containers that have broken in transit. Specimens received beyond acceptable holding times (see Specimen Requirements). Specimens collected after the expiration date on the collection kit. Specimens collected in collection devices from other manufacturers. Incorrect volume of urine in urine transport tube.
<b>Requisition Form</b>	Clinical Test Requisition (select <b>Chlamydia &amp; Gonorrhea Nucleic Acid Amplification Test</b> ).
<b>Required Information</b>	Name and address of submitter (and/or Horizon profile #). Patient name or identifier, town of residence (city, state, zip), date of birth. Specimen source/type, date collected and test requested.
<b>Limitations</b>	The test is not intended for the evaluation of suspected sexual abuse or for other medico-legal indications. Therapeutic failure or success cannot be determined with this test since nucleic acid may persist following appropriate antimicrobial therapy. Test results should be interpreted in conjunction with all other laboratory and clinical data available to the clinician.