

<b>CDC Influenza SARS-CoV-2 Multiplex (Emergency Use Authorization)</b>	
<b>Test Description</b>	Assay for the simultaneous qualitative detection and differentiation of nucleic acids from the SARS-CoV-2, Influenza A, and/or Influenza B viruses in upper or lower respiratory specimens
<b>Test Use</b>	To aid in the diagnosis of individuals suspected of respiratory viral infection consistent with COVID-19 for public health surveillance purposes or outbreak investigations
<b>Test Department</b>	Advanced Molecular Diagnostics Phone: (860) 920-6662, FAX: (860) 920-6661
<b>Methodology</b>	Real-time RT-PCR multiplexed test
<b>Availability</b>	Daily, Monday-Friday
<b>Specimen Requirements</b>	Nasopharyngeal, oropharyngeal, nasal, or mid-turbinate swab specimen collected and placed into a transport tube containing 2-3 mL viral or universal transport medium, Amies transport medium (eSwab), or sterile saline. Sputum, BAL or nasopharyngeal aspirate, and the non-bacteriostatic saline used to collect the specimen should be placed immediately into a sterile transport tube
<b>Collection Kit/Container</b>	Specimen collection kits can be obtained by calling the SPHL outfit room, (860) 920-6674 or (860) 920-6675, Monday- Friday, 8:00 AM to 4:00 PM. Requests may also be submitted via e-mail to <a href="mailto:dph.outfitroom@ct.gov">dph.outfitroom@ct.gov</a>
<b>Collection Instructions</b>	Use freshly collected specimens for optimal test performance. Inadequate specimen collection or improper sample handling, storage, or transport may yield erroneous results. Refer to the CDC Interim Guidelines for Collecting, Handling, and Testing Clinical Specimens from Persons for Coronavirus Disease 2019 (COVID-19) <a href="https://www.cdc.gov/coronavirus/2019-nCoV/lab/guidelines-clinical-specimens.html">https://www.cdc.gov/coronavirus/2019-nCoV/lab/guidelines-clinical-specimens.html</a>
<b>Specimen Handling &amp; Transport</b>	Store specimen at 2-8°C up to 3 days. Transport to the laboratory with a frozen ice pack coolant. If there is a delay in shipment expected, store specimens at -70°C or lower until delivered to the laboratory.
<b>Unacceptable Conditions</b>	Unlabeled specimen; Improper specimen type; Specimens that have leaked or containers that have broken in transit; Specimens not handled, stored, or transported as described above.
<b>Requisition Form</b>	Indicate <b>Flu/SARS-CoV-2 or Flu/COVID-19 PCR</b> in the "Test, Agent or Disease Not Listed" box of the OL9B Clinical Test Requisition form.
<b>Required Information</b>	Name and address of submitter (and/or Horizon profile #) Patient name or identifier, town of residence (city, state, zip), date of birth Specimen type or source of collection, date collected, test requested Please ensure patient name on the requisition matches that on the specimen.
<b>Limitations</b>	Performance of the Influenza SARS-CoV-2 (Flu SC2) Multiplex Assay has only been established in upper and lower respiratory specimens. Negative results do not preclude SARS-CoV-2 infection and should not be used as the sole basis for treatment or other patient management decisions. Performance evaluations have shown decreased sensitivity of the influenza A target when a high titer of SARS-CoV2 or influenza B is also present in the sample.
<b>Additional comments</b>	The possibility of a false negative result should especially be considered if the patient's recent exposures or clinical presentation suggest that SARS-CoV-2 infection is possible. In the case of influenza A and B viruses, children tend to shed virus more abundantly and for longer periods of time than adults. The performance of the assay has not been established in individuals who received nasally administered influenza vaccine.

Revision: 09/28/2020