

Document Title: Screening Test for Semen

Controlled: Yes, with red stamp present

Controlled By: Quality Manager

Prepared By: _____ Date: _____

Approved By: _____ Date: _____

A. PURPOSE:

To perform a screening test for the presence of semen in Forensic samples.

B. RESPONSIBILITY:

Forensic Science Examiners from the Connecticut State Forensic Science Laboratory who have been trained in the discipline of semen screening according to SOP-FB-31 (Training Manual).

C. DEFINITIONS:

1. PBS: Phosphate Buffered Saline
2. AP: Acid Phosphatase

D. PROCEDURE:

This test will be performed at the discretion of the examiner based on the submitting agency requests, case information and the condition of the evidence.

1. Materials:
 - a. Acid Phosphatase Reagent
 - aa. α -naphthyl phosphate substrate solution
 - bb. Fast Blue B color solution
 - b. Controls: positive (known 1:10 semen stain) and negative (blank filter paper)
 - c. PBS
 - d. Cotton swabs or spot plates
2. Procedure:
 - a. Test a positive control and a negative control (blank filter paper) with PBS according to the following procedure (steps 2.b. – 2.e.) prior to the questioned sample.
 - aa. If controls yield the appropriate results, record on the appropriate Quality Record Worksheet and test the questioned samples.
 - bb. If controls do not yield the appropriate results, review the procedure and retest the controls prior to the questioned samples.
 - cc. If necessary, the reagent may be retested with the controls within the workday.

- D. 2. b. Remove and test a portion of the questioned stain or area.
 - c. Add one drop of α -naphthyl phosphate substrate reagent to the sample. Wait 10 seconds.
 - d. Add one drop of Fast Blue B color reagent to the α -naphthyl phosphate reagent and sample.
 - e. Observe any color change within 15 seconds.
 - f. Discard any unused reagent daily.
3. Results:
 - a. *Positive*. The development of a purple or pink color within 15 seconds indicates a positive result and detects the presence of acid phosphatase activity.
 - b. *Negative*. No color change indicates a negative result and no acid phosphatase activity is detected.
 - c. *Inconclusive*. No discernible color change.
 - d. Record the results on the appropriate Quality Record Worksheet.
 4. Record reagent used on the General Reagent Sheet (FBQR-09).

E. REFERENCES:

1. Blake, E.T. and G.F. Sensabaugh: "Genetic markers In Human Semen: A Review". Journal of Forensic Sciences, Vol. 21, 784-796, 1976.
2. Gaensslen, R. E. , Sourcebook In Forensic Serology, Immunology, and Biochemistry , U.S. Government Printing Office, Washington D.C., 1983.
3. Metropolitan Police Forensic Science Laboratory. Biology Methods Manual. 1978, pp.3-17 to 3-20.