

DNA Training Manual

Individualized training plan prepared for: _____

Date of start of training: _____

Date of completion of bench work training: _____

Date of completion of analysis and report writing training: _____

Date of Technical Review authorization: _____

7.1 Purpose

The primary objective reflected in this section is to demonstrate a comprehensive training program for all personnel entering the Division of Scientific Services DNA Unit, or being trained in particular tasks in the DNA Unit, to allow for, at its completion, employees to be fully capable of independent casework dependent upon their job title and assigned duties.

This program is designed to train personnel in laboratory-specific forensic DNA protocols. Here within is an outline of formal procedures for the training and assessment of new examiners and technicians in the DNA unit. The goal of this program is to develop DNA examiners and/or technicians capable of performing independent forensic DNA analysis. Successful completion of this program in its entirety will take approximately six months. However, the training program may be abbreviated as warranted for examiners and technicians that have previous experience in forensic DNA methods with approval of the TL. Each trainee will be given this training manual with the tasks that they are to complete documented by a check in the box to the left of each task. While new employees just starting their careers in forensic DNA examination will need to complete all tasks in their assigned modules, new employees with previous forensic DNA experience will have modified training plans geared due to their previous expertise. Therefore, some listed tasks will not be required, and will therefore not be checked off. Also, not every employee will be trained in every module. Some new employees will not be tasked with, for example, database processing or known DNA analysis. Therefore, it is possible that full modules of this manual will remain blank. Regardless of full or modified training, a competency test will be given prior to performing casework. The signature of the technical leader, assistant director, and deputy director on the last page of this manual will be documentation of his/her approval of the individual training plan, including the approval of omitted parts of the plan due to analysts previous experience.

This training program is designed to supplement, for analysts, successful completion of college coursework in biochemistry, molecular biology, genetics, and population genetics/statistics as required by the FBI QAS. The training program will expose the DNA trainee to the scientific principles underlying each DNA test used by the State of Connecticut Division of Scientific Services. Prior to the administering of competency tests and moot courts, a review of all training documents is to be completed by the Technical Leader, Quality Manager, and Assistant Director and/or Deputy Director. Moot courts will include direct and cross examination, as well as the introduction of evidence/exhibits. There will be documentation of the moot court exercise, containing an evaluation of the analyst's performance that will be retained by the laboratory. Assessment will be made at the completion of each module. Each module will be approved by the Training Coordinator, supervisor, or Training Coordinator designee once the trainee has been deemed competent at the specified tasks. Multiple modules may be assessed simultaneously.

The trainee will maintain (a) training binder(s) containing all worksheets, printouts generated, and results of each module as they are being completed. At the completion of training, the Training Coordinator will maintain the binder(s) in his/her office, being retained at the laboratory (please refer to GL-11 for the control of records for record retention policy). These records should clearly indicate which training module they pertain to. Most DNA Unit personnel must pursue yearly continuing education in the form of workshops, seminars, professional conferences or collegiate coursework as well as further on-the-job training as specified by FBI QAS, to be documented separate from this training binder.

If additional training above those tasks listed in this manual becomes necessary, a modified training plan must be created, documented on DNA QR-283 Addendum Training Record, and approved by the Technical Leader prior to implementation. A competency test will still be administered at the completion of this addendum.

For any necessary re-training of DNA Unit personnel, please refer to GL-14 General Training.

7.2 Responsibility

DNA Unit personnel: New hires to complete, and current personnel to assist in/augment their training process.

7.3 DNA Training Program: New Employees or Additional Training, DNA Unit

7.3.1 Introduction

This section must be completed by Forensic Science Examiners 1, 2, and 3, Laboratory Assistants, and Connecticut Careers Trainees (CCT).

7.3.1.1 Goal: Upon completion, the examiner will be familiar with the general forensic laboratory operations and his/her individual responsibilities.

7.3.1.2 Tasks

7.3.1.2.1 ☐ Orientation to the laboratory facility, personnel, table of organization, and the chain of command

7.3.1.2.2 ☐ Familiarization with other forensic disciplines practiced at the laboratory, via laboratory tour

Trainee Initials/Date	Trainer Initials/Date
--------------------------	--------------------------

--	--

--	--

Approved by Director: Dr. Guy Vallaro

7.3.1.3	Required Reading	Trainee Initial/Date	Trainer Initial/Date
7.3.1.3.1	Laboratory Quality Manual <input type="checkbox"/>		
7.3.1.3.2	SWGDM Guidelines – current version <input type="checkbox"/>		
7.3.1.3.3	Quality Assurance Standards for Forensic DNA Testing Laboratories (current version) <input type="checkbox"/>		
7.3.1.3.4	Quality Assurance Standards for DNA Databasing Laboratories (current version) <input type="checkbox"/>		
7.3.1.4	Assessment of Introduction Section		
7.3.1.4.1	Oral and/or written evaluation by the supervisor or designee, to include moot court questions <input type="checkbox"/>		
7.3.2	Evidence Handling / Examination This section must be completed by Forensic Science Examiners 1, 2, and 3, Laboratory Assistants, and CCTs.		
7.3.2.1	Goals		
7.3.2.1.1	To handle evidentiary samples in an appropriate manner		
7.3.2.1.2	To preserve evidence that may need to be analyzed by other sections within the laboratory		
7.3.2.1.3	To learn the operation of the LIMS computer system, understanding the importance of the chain of custody, and the creation of items of evidence at the time they are made.		
7.3.2.1.4	To demonstrate competency in the basic tasks necessary to complete evidence documentation and handling for DNA (DNA QR-1 and DNA-QR-1A) through written exam and/or questions given in a moot-court setting.		
7.3.2.2	Tasks		

DNA SOP-7 DNA Training Manual*Approved by Director: Dr. Guy Vallaro*

Document ID: 927

Revision: 18

Effective Date: 8/20/2021

Status: Published

Page 5 of 34

		Trainee Initials/Date	Trainer Initials/Date
7.3.2.2.1	Receive training on LIMS computer system from Laboratory LIMS Administrator		
<input type="checkbox"/>			
7.3.2.2.2	To learn the laboratory procedures for receipt, transfer, storage, and return of evidence on the LIMS computer system for DNA casework		
<input type="checkbox"/>			
7.3.2.2.3	To learn the written procedures and LIMS computer system for receipt, transfer, and storage of convicted offender samples for DNA database analysis		
<input type="checkbox"/>			
7.3.2.2.4	To learn how and when to create sub-items for evidence on the LIMS computer system		
<input type="checkbox"/>			
7.3.2.2.5	To demonstrate knowledge of how evidence is stored at the laboratory		
<input type="checkbox"/>			
7.3.2.2.6	To demonstrate knowledge of safe handling procedures of evidence (to avoid contamination of evidence or exposure of examiner/co-workers to potential biohazards)		
<input type="checkbox"/>			
7.3.2.2.7	To demonstrate knowledge of maintaining the the custody and integrity of evidence.		
<input type="checkbox"/>			
7.3.2.2.8	To learn the laboratory case acceptance policies by reading ER SOP-01, ER SOP 1.1, ER SOP-02, ER SOP 02-1, and ER SOP-04		
<input type="checkbox"/>			
7.3.2.2.9	To demonstrate knowledge of awareness of order of evidence examinations.		
<input type="checkbox"/>			
7.3.2.2.10	Orientation to Evidence Receiving, given by an ECO		
<input type="checkbox"/>			

Approved by Director: Dr. Guy Vallaro

7.3.2.3 Required Reading

Trainee Initials/Date	Trainer Initials/Date

7.3.2.3.1 Quality Manual for LIMS (GL-4)

☐

7.3.2.4 Assessment of Evidence Handling section

7.3.2.4.1 Oral and/or written evaluation by the supervisor

☐

Or designee, to include moot court questions

7.3.3 Foundational Scientific Knowledge

This section must be completed by Forensic Science Examiners 1, 2, and 3, and CCTs. Laboratory assistants may complete tasks deemed appropriate by the supervisor, trainer, or designee.

7.3.3.1 Goal: To ensure examiners have both the formal education and a working knowledge of the fundamental scientific concepts underlying DNA analysis in a forensic laboratory setting.

7.3.3.2 Tasks: To document coursework and/or training in the following general subject areas and a working knowledge of the principles listed, providing transcripts. This approval shall be completed by the Technical Leader after review of transcripts.

7.3.3.2.1 Genetics: Alleles; Mendelian inheritance; genotype vs. phenotype; coding vs. non-coding, DNA vs. protein markers

☐

Trainee Initials/Date	Technical Leader Initial/Date

7.3.3.2.2 Biochemistry: The study of the nature of biologically important molecules in living systems, DNA replication and protein synthesis, and the quantitative and qualitative aspects of cellular metabolism

☐

7.3.3.2.3 Molecular Biology: the study of theories, methods, and techniques used in the study and analysis of gene structure, organization, and function

☐

Approved by Director: Dr. Guy Vallaro

		Trainee Initials/Date	Technical Leader Initials/Date
7.3.3.2.4	Population genetics and statistics: Hardy-Weinberg equilibrium; allelic and genotypic frequencies; ideal population considerations; theta, population databases and minimum size; the use of different race/ethnicity population databases; population substructure; expected frequency (match probability) vs. likelihood ratio calculations; basic probability concepts (product rule, independence)		
<input type="checkbox"/>			
7.3.3.2.5	Document in training binder relevant continuing education and/or training (in-house or outside agency), providing certificates and/or topic outlines		
<input type="checkbox"/>			
7.3.3.3	Required Reading		
7.3.3.3.1	NRC II (1996)		
<input type="checkbox"/>			
7.3.3.3.2	Advanced Topics in Forensic DNA Typing: Methodology (Copyright 2012). (Selected chapters, as determined by trainer)		
<input type="checkbox"/>			
7.3.3.4	Assessment of Foundational Scientific Knowledge Section		
7.3.3.4.1	Oral examination in moot court setting covering basic principles and required reading		
<input type="checkbox"/>			

7.3.4 Applied Scientific Knowledge

This section must be completed by Forensic Science Examiners 1, 2, and 3 and CCTs. Laboratory assistants may complete tasks deemed appropriate by the supervisor or designee.

Approved by Director: Dr. Guy Vallaro

7.3.4.1 Goal: To train an examiner in the details of forensic DNA analyses such that he/she will be able to (a) Apply the knowledge to the processing of forensic DNA evidence and (b) Provide the court with an appropriate explanation of how DNA testing is performed.

7.3.4.2 Tasks: to provide instruction and theory regarding the following:

		Trainee Initials/Date	Trainer Initials/Date
7.3.4.2.1	Basic Biochemical Formulas:		
<input type="checkbox"/>	Performing serial dilutions, determining concentrations, etc.		
7.3.4.2.2	Basic theoretical understanding of past methodologies and platforms of forensic DNA typing (through assigned chapter readings)		
<input type="checkbox"/>			
7.3.4.2.3	DNA Extraction Methods: Provide knowledge of the protocol differences for DNA extractions (manual and automation) of blood, buccal cells, hair, bone, teeth, tissue, and differential DNA extractions		
<input type="checkbox"/>			
7.3.4.2.4	DNA Quantitation Method: Quantifiler Trio using 7500 real-time PCR and software		
<input type="checkbox"/>			
7.3.4.2.5	PCR Based Methods		
7.3.4.2.5.1	Autosomal STRs		
<input type="checkbox"/>			
7.3.4.2.5.2	Y-STRs		
<input type="checkbox"/>			
7.3.4.2.6	Population Statistics		
7.3.4.2.6.1	Determining Allele Frequencies		
<input type="checkbox"/>			
7.3.4.2.6.2	National Population DNA Databases		
<input type="checkbox"/>			

		Trainee Initials/Date	Trainer Initials/Date
7.3.4.3	Required Reading		
7.3.4.3.1 <input type="checkbox"/>	Internal validation summaries of currently used methodologies performed by the Division of Scientific Services		
7.3.4.3.2 <input type="checkbox"/>	Developmental validation studies of amplification kits currently used by the Division of Scientific Services		
7.3.4.3.3 <input type="checkbox"/>	Review articles for STRs (Autosomal & Ys)		
7.3.4.3.2 <input type="checkbox"/>	SWGDM Contamination Prevention and Detection Guidelines for Forensic DNA Laboratories		
7.3.4.3.3 <input type="checkbox"/>	Familiarization with risk assessment (noted in GL-1 and GL-9, briefly, and prevention of sample loss through reducing/preventing contamination		
7.3.4.4	Assessment of Applied Scientific Knowledge section		
7.3.4.4.1 <input type="checkbox"/>	Moot court to demonstrate theoretical understanding of all tasks/readings in section		

7.3.5 Laboratory Analytical Procedures

Definitions:

Trainer: Laboratory personnel competent in specific workflow process being taught.

Trainee: Laboratory personnel learning the workflow process.

Observed: Trainee will observe the trainer perform the workflow process.

Supervised: Trainer will observe the trainee perform the workflow process.

Independent: Trainee will perform the workflow process without supervision.

7.3.5.1 Casework Laboratory Procedures

This section must be completed by Forensic Science Examiners 1, 2, and 3, laboratory assistants, and CCTs working with casework evidentiary samples.

Extracted samples will be quantitated, amplified, run on a genetic analyzer, and analyzed to determine if a DNA profile is detectable from the extracted genomic DNA.

7.3.5.1.1 Goal: To provide practical instruction on the laboratory procedures to be used by the examiner, ending with a competency test covering all aspects and procedures trained in.

Note: Completion of Training: The trainer and trainee can evaluate the training and determine if additional or less training in a specific area is necessary. This alteration to the training must be documented and agreed upon by training coordinator, trainer, trainee, and Technical Leader. Examination and competency test will not be taken until training coordinator, trainer, and trainee deem the trainee ready.

7.3.5.1.2	Tasks	Trainee Initials/Date	Trainer Initials/Date
7.3.5.1.2.1 <input type="checkbox"/>	SOP Review: DNA SOP-20 Extraction of Unknown Samples on EZ1 Advanced XL		
7.3.5.1.2.2 <input type="checkbox"/>	Review Reagent QC List		
7.3.5.1.2.3 <input type="checkbox"/>	Training in DNA Worklist macro, to make worksheets for extraction and quantification		
7.3.5.1.2.4 <input type="checkbox"/>	(To be conducted by a member of Forensic Biology section) KM training, including observed and supervised testing of various stains and substrates, and theory.		
7.3.5.1.2.5 <input type="checkbox"/>	(To be conducted by a member of Forensic Biology section) Evidence examination training of swab-only evidence submissions		

Approved by Director: Dr. Guy Vallaro

		Trainee Initials/Date	Trainer Initials/Date
<input type="checkbox"/>	7.3.5.1.2.6 Understand the difference between evidence and work product		
<input type="checkbox"/>	7.3.5.1.2.7 Reading: Product inserts/instruction manuals for EZ1 and Investigator Kits		
<input type="checkbox"/>	7.3.5.1.2.8 Observed lab work: DNA extraction of real or mock evidence on Qiagen EZ1: at least three samples, including one hair, appropriate positive controls and reagent blanks		
<input type="checkbox"/>	7.3.5.1.2.9 Review PowerPoint presentation on theory of qPCR		
<input type="checkbox"/>	7.3.5.1.2.10 SOP Review: DNA SOP-3 DNA Quantitation		
<input type="checkbox"/>	7.3.5.1.2.11 Reading: product insert/instruction manuals for 7500 and Quantifiler Trio		
<input type="checkbox"/>	7.3.5.1.2.12 SOP Review: DNA WI-07 Quantifiler Trio DNA Quantification Work Instructions, paying attention to "stop at quant" procedure, and how degradation affects amplification.		
<input type="checkbox"/>	7.3.5.1.2.13 Training in 7500 software setup, to include use of macros		
<input type="checkbox"/>	7.3.5.1.2.14 Observed lab work: DNA quantitation of samples previously extracted		
<input type="checkbox"/>	7.3.5.1.2.15 Software Training: 7500 software, analysis, & troubleshooting		

Approved by Director: Dr. Guy Vallaro

		Trainer Initials/Date	Trainee Initials/date
7.3.5.1.2.16	Training in amplification, dilution, halt at quant, and staff search, worksheet creation		
<input type="checkbox"/>			
7.3.5.1.2.17	SOP Review: SOP for current STR kit that covers its amplification and detection		
<input type="checkbox"/>			
2.3.5.1.2.18	Overview of history of legacy typing kits		
<input type="checkbox"/>			
7.3.5.1.2.19	Reading: Product insert for current STR amplification system		
<input type="checkbox"/>			
7.3.5.1.2.20	Reading: Product insert for current Y-STR amplification system		
<input type="checkbox"/>			
7.3.5.1.2.21	Observed lab work: Amplification of samples with STR kit currently used in casework		
<input type="checkbox"/>			
7.3.5.1.2.22	Observed lab work: Amplification of appropriate samples with Y-STR kit currently used for casework		
<input type="checkbox"/>			
7.3.5.1.2.23	Workbook macro training to create CE injection plate		
<input type="checkbox"/>			
7.3.5.1.2.24	Learn routine maintenance of CE instrument		
<input type="checkbox"/>			
7.3.5.1.2.25	Observed lab work: Operation of CE collection software, preparation of an injection plate for samples previously amplified		
<input type="checkbox"/>			
7.3.5.1.2.26	Supervised lab work: DNA Extraction, quant, amplification (STR & Y-STR) through injection of mock evidence-type samples: five samples, one positive control, and one reagent blank. Creation of worksheets need not be supervised		
<input type="checkbox"/>			

Approved by Director: Dr. Guy Vallaro

		Trainee Initials/Date	Trainer Initials/Date
7.3.5.1.2.27	Independent lab work: DNA Extraction, quant, amplification (STR& Y-STR) through injection of mock evidence-type samples: ten samples (including two hairs) , one positive control, and one reagent blank, including creation of worksheets, and importing samples into appropriate analysis software. (2X)	1) <input type="text"/>	<input type="text"/>
<input type="checkbox"/>		2) <input type="text"/>	<input type="text"/>
7.3.5.1.2.28	Observed lab work: DNA Extraction of semen-containing evidence samples on Qiagen EZ1: one Sample, one positive control, and one reagent blank. (If previous EZ1 run has been observed, only necessary for observation through incubation of b-fractions.)	<input type="text"/>	<input type="text"/>
<input type="checkbox"/>			
7.3.5.1.2.29	Supervised lab work: DNA Extraction of semen-containing evidence samples on Qiagen EZ1: two samples, one positive control, and one reagent blank, through injection. (If previous EZ1 run has been supervised, only necessary to supervise through incubation of b-fractions.)	<input type="text"/>	<input type="text"/>
<input type="checkbox"/>			
7.3.5.1.2.30	Independent lab work: DNA Extraction of semen-containing evidence samples on Qiagen EZ1: five samples one positive control, and one reagent blank, through injection. (2X)	1) <input type="text"/>	<input type="text"/>
<input type="checkbox"/>		2) <input type="text"/>	<input type="text"/>
7.3.5.1.2.31	SOP Review: DNA SOP-20.6 Concentration and Purification of DNA Samples, to include discussion on how to determine Microcon appropriateness.	<input type="text"/>	<input type="text"/>
<input type="checkbox"/>			
7.3.5.1.2.32	Observed lab work: Concentration of an appropriate sample previously extracted, or of a mock sample (will not proceed past concentration).	<input type="text"/>	<input type="text"/>
<input type="checkbox"/>			

Approved by Director: Dr. Guy Vallaro

		Trainee Initials/Date	Trainer Initials/Date
7.3.5.1.2.33	Supervised lab work: Concentrate an appropriate sample previously extracted		
<input type="checkbox"/>			
7.3.5.1.2.34	Independent lab work: Concentrate 2 previously extracted samples		
<input type="checkbox"/>			
7.3.5.1.2.35	Independent lab work: Perform amplification through injection of 3 concentrated samples		
<input type="checkbox"/>			
7.3.5.1.2.36	Male Screen Training		
7.3.5.1.2.36.1	Review of QC, worksheets, & LIMS requests necessary to perform procedure		
<input type="checkbox"/>			
7.3.5.1.2.36.2	Review of appropriate validation materials		
<input type="checkbox"/>			
7.3.5.1.2.36.3	Observed work: Demonstration of male screen extraction procedure of four swabs, each containing 90-140ul female saliva (depending on swab size) and 10ul semen (neat, 1:10 dil, 1:100 dil & 1:1000 dil). Trainer selected samples will be Microcon concentrated.		
<input type="checkbox"/>			
7.3.5.1.2.36.4	Supervised work: Male screen extraction procedure of duplicated four swabs above. Trainer selected samples will be Microcon concentrated. Samples from observed & supervised extractions quanted by trainer to show similarity.		
<input type="checkbox"/>			
7.3.5.1.2.36.5	Independent work: male screen extraction procedure of two additional sets identical to sets previously processed. Samples quanted in duplicate. Trainer selected samples will be Microcon concentrated. Results compared to two previous sets. Sets processed simultaneously, each set including appropriate controls.		
<input type="checkbox"/>			

Approved by Director: Dr. Guy Vallaro

		Trainee Initials/Date	Trainer Initials/Date
7.3.5.1.2.36.6	Competency Test: Male Screen Procedure, to include practical, written and oral assessments		
<input type="checkbox"/>			
7.3.5.1.3	Assessment of Section Casework Laboratory Procedures section: Competency Exam to include:		
7.3.5.1.3.1	<u>Evidentiary Sample Practical</u> : To include evidence six non-semen containing (including two hairs), and two semen-containing samples, with appropriate controls, EZ1 extraction through injection. A minimum of one sample to be concentrated, and a minimum of one sample to be amplified with Y-STRs		
<input type="checkbox"/>			
7.3.5.1.3.2	<u>Written Exam</u> : pertaining to evidentiary sample processing		
<input type="checkbox"/>			
7.3.5.1.3.3	<u>Moot court</u> : covering laboratory techniques from evidence examination through injection		
<input type="checkbox"/>			

7.3.5.2 Database/Casework Knowns Laboratory Procedures

This section must be completed by Forensic Science Examiners 1, 2, and 3, Laboratory Assistants, and CCTs working with known evidentiary and/or database samples.

Extracted samples will be amplified, injected, and analyzed to determine if correct DNA profile is generated from the extracted genomic DNA. Convicted offender training samples should all be in CODIS with expanded loci. Competency test samples will not be samples used in training.

- 7.3.5.2.1 Goal: To provide practical instruction on the laboratory procedures to be used by the examiner concluded with a competency test, to cover all aspects and procedures in which one is trained.

Note: Completion of Training: The trainer and trainee can evaluate the training and determine if additional or less training in a specific area is necessary. This alteration to the training must be documented and agreed upon by training coordinator, trainer, trainee, and Technical Leader. Examination and competency test will not be taken until training coordinator, trainer, and trainee deem the trainee ready.

Approved by Director: Dr. Guy Vallaro

7.3.5.2.2 Tasks		Trainee Initials/Date	Trainer Initials/Date
7.3.5.2.2.1 <input type="checkbox"/>	SOP Review: DNA SOP-19 Processing of Single Source Samples on the EZ1 Advanced XL		
7.3.5.2.2.2 <input type="checkbox"/>	Reading: Product inserts/instruction manuals for EZ1 and Investigator Kits		
7.3.5.2.2.3 <input type="checkbox"/>	Observed lab work: Extraction of one blood, one buccal FTA, one buccal swab, one positive control and one reagent blank using normalization procedure on EZ1(to include worksheet creation)		
7.3.5.2.2.4 <input type="checkbox"/>	Reading: Product inserts/instruction manuals for STR & Y-STR amplification kits		
7.3.5.2.2.5 <input type="checkbox"/>	Observed lab work: Amplification (STRs & Y-STRs) of samples previously extracted (to include worksheet creation)		
7.3.5.2.2.6 <input type="checkbox"/>	Learn routine maintenance of CE instrument		
7.3.5.2.2.7 <input type="checkbox"/>	Observed lab work: Operation of CE instrument collection software, preparation of an injection plate for samples previously amplified (to include worksheet creation)		
7.3.5.2.2.8 <input type="checkbox"/>	Supervised lab work: Extraction of 3 bloods, 2 Buccal FTAs, 1 buccal swab, and appropriate controls using normalization procedure on EZ1, manual amplification (STRs & Y-STRs) and injection		
7.3.5.2.2.9 <input type="checkbox"/>	Independent lab work: Extraction of 5 bloods, 3 Buccal FTAs, 2 buccal swabs, and appropriate controls using normalization procedure on EZ1, manual amplification (STRs and Y-STRs) and injection		

Approved by Director: Dr. Guy Vallaro

		Trainee Initials/Date	Trainer Initials/Date
7.3.5.2.2.10	SOP review: SOPs pertaining to processing of casework knowns and single source samples		
<input type="checkbox"/>			
7.3.5.2.2.11	Use macro to create sample list of consecutive database or known samples, and setup plate workbook		
<input type="checkbox"/>			
7.3.5.2.2.12	Reading: BSD Duet users' manual		
<input type="checkbox"/>			
7.3.5.2.2.13	Training on BSD cleaning & Maintenance		
<input type="checkbox"/>			
7.3.5.2.2.14	Observed lab work: Set-up & punching of up to one full plate		
<input type="checkbox"/>			
7.3.5.2.2.15	Reading: GlobalFiler Express insert		
<input type="checkbox"/>			
7.3.5.2.2.16	Observed lab work: Direct Amplification of up to one full plate of samples previously punched		
<input type="checkbox"/>			
7.3.5.2.2.17	Observed lab work: Injection of up to one full plate of samples, previously amplified		
<input type="checkbox"/>			
7.3.5.2.2.18	Observed lab work: Reinjections (cherry-picking) of samples previously injected (minimum of 5 samples)		
<input type="checkbox"/>			
7.3.5.2.2.19	Supervised lab work: Setup & punching of plate containing at least two full modules, followed by direct amplification, injection, and cherry-picking of samples needing reinjection (determined by analyst, minimum of 5)		
<input type="checkbox"/>			

Approved by Director: Dr. Guy Vallaro

		Trainee Initials/Date	Trainer Initials/Date
7.3.5.2.2.20	Independent lab work: Setup & punching of plate containing at least two full modules, followed by direct amplification, injection, and cherry-picking of samples needing reinjection (determined by analyst, minimum of 5) (X2)	1)	
<input type="checkbox"/>		2)	
7.3.5.2.2.21	Macro training to create sample list of non-consecutive database or known samples, and setup plate workbook		
<input type="checkbox"/>			
7.3.5.2.2.22	Observed lab work: Entire process of known examination, from evidence receiving, to storage, to examiner, to exam, back to storage. To include write-up and LIMS itemization. At least 10 knowns, including 2 FTA cards, 2 swabs, one blood FTA and one EZ collect device.		
<input type="checkbox"/>			
7.3.5.2.2.23	Supervised lab work: Evidence examination of 2 mock knowns (one FTA and one buccal swab)		
<input type="checkbox"/>			
7.3.5.2.2.24	Independent lab work: Evidence examination of 2 mock knowns (one FTA and one buccal swab)		
<input type="checkbox"/>			
7.3.5.2.2.25	Competency Exam: Set-up through injection of one module, using direct amplification, and 2 blood & 2 buccal samples on the EZ1 using the normalization procedure with moot court to demonstrate understanding of laboratory processes & concepts, as well as a written exam		
<input type="checkbox"/>			

7.3.5.3 Convicted Offender Samples and CODIS

This section must be completed by Forensic Science Examiners 1, 2, and 3, Laboratory Assistants and CCTs working with convicted offender samples. It should be completed in addition to QR-284, CODIS Training.

Approved by Director: Dr. Guy Vallaro

		Trainee Initials/Date	Trainer Initials/Date
7.3.5.3.1 <input type="checkbox"/>	SOP Review: DNA SOP-10, CODIS Administration, DNA SOP-15, CODIS Hit Confirmations, DNA SOP-11, Collection from Offenders, Missing Persons and Relatives of Missing		
7.3.5.3.2 <input type="checkbox"/>	Reading: State database statutes (CODIS)		
7.3.5.3.3 <input type="checkbox"/>	Reading: State database regulations (CODIS)		
7.3.5.3.4 <input type="checkbox"/>	CODIS Training on CJIS-WAN (modules taken at the discretion of the CODIS administrator)		
7.3.5.3.5 <input type="checkbox"/>	Observation: import of one full plate and one partial late of database samples including checking on duplicates & dispositions		
7.3.5.3.6 <input type="checkbox"/>	Supervised work: import of one full plate and one partial plate of database samples, including checking on duplicates & dispositions.		
7.3.5.3.7 <input type="checkbox"/>	Independent work: import of one full plate and one partial plate of database samples, including checking on duplicates & dispositions.		
7.3.5.3.8 <input type="checkbox"/>	COLLECT training		
7.3.5.3.9 <input type="checkbox"/>	Training on hit confirmations: Understanding workflow to check for new hit confirmations, pull database cards, COLLECT searches, and processing of hit confirmation paperwork		

Approved by Director: Dr. Guy Vallaro

		Trainee Initials/Date	Trainer Initials/Date
7.3.5.3.10	Training on post-processing of convicted offender samples		
<input type="checkbox"/>			
7.3.5.3.11	Observed work: post-processing of 10 samples		
<input type="checkbox"/>			
7.3.5.3.12	Supervised work: post-processing of 20 samples		
<input type="checkbox"/>			
7.3.5.3.13	Competency test: Importing of one full plate of database samples, written exam, & moot court Describing laboratory CODIS policies.		
<input type="checkbox"/>			
7.3.5.3.14	DNA QR-284, CODIS Training, completed		
<input type="checkbox"/>			
7.3.5.4	Analysis for Questioned & Known Samples		
	To be completed by Forensic Science Examiners 1, 2, and 3 and CCTs (prior to completion of program), who will be doing complex analysis on all types of forensic samples, including STRmix interpretation, with the eventual goal to be able to seamlessly portray all analyses and protocols in a court of law to a lay jury member. Sections pertaining to administrative review can be completed without the need for analysis or technical review training. Employees who are not deemed qualified analysts may be trained in, and complete administrative reviews of both batch paperwork and case jackets.		
7.3.5.4.1	SOP Review: DNA SOPs pertaining to analysis of STRs and Y-STRs, including software programs.		
<input type="checkbox"/>			
7.3.5.4.2	GeneMarker Demonstration		
<input type="checkbox"/>			
7.3.5.4.3	Review PowerPoint presentation on common artifacts found in forensic analysis		
<input type="checkbox"/>			

Approved by Director: Dr. Guy Vallaro

		Trainee Initials/Date	Trainer Initials/Date
7.3.5.4.4	Reading: SWGDAM Interpretation Guidelines for Autosomal STR Typing by Forensic DNA Testing Laboratories		
<input type="checkbox"/>			
7.3.5.4.5	Number of contributor training, including training on theory, utilization, and limitations of DNA-QR-302 Contributor Estimation Worksheet		
<input type="checkbox"/>			
7.3.5.4.6	Observed analysis: Review of a batch through analysis (3X) (to include macros used during analysis process such as concordance check, project comparison, etc.)	1)	
<input type="checkbox"/>		2)	
		3)	
7.3.5.4.7	Review case jacket documentation upon batch completion		
<input type="checkbox"/>			
7.3.5.4.8	Observed analysis: Technical Review of a batch through analysis (3X)	1)	
<input type="checkbox"/>		2)	
		3)	
7.3.5.4.9	Observed analysis: Administrative Review of a batch through analysis		
<input type="checkbox"/>			
7.3.5.4.10	Independent analysis: review of mock-batch (to include a minimum of 30 STR profiles) for analysis, technical review, or administrative, review, as appropriate for required duties (3X)	1)	
<input type="checkbox"/>		2)	
		3)	
7.3.5.4.11	Demonstration on GeneMapper IDX (Y-STR analysis only)		
<input type="checkbox"/>			
7.3.5.4.12	Reading: SWGDAM Interpretation Guidelines for Y-Chromosome STR Typing		
<input type="checkbox"/>			

Approved by Director: Dr. Guy Vallaro

Trainee Initials/Date	Trainer Initials/Date
--------------------------	--------------------------

7.3.5.4.13 Observed analysis: 5 Y-STR profiles

☐

7.3.5.4.14 Independent analysis: analysis of Y-STR data (at least 5 questioned samples)

☐

7.3.5.4.15 SOP Review: SOPs pertaining to processing of convicted offender/database/single source samples using current methods

☐

7.3.5.4.16 GeneMarker demonstration, specifically for known and/or database sample analysis

☐

7.3.5.4.17 Macro training for all macros associated with known analysis, including staff search when contamination is detected

☐

7.3.5.4.18 Observed analysis: analysis of one batch of casework known samples

☐

7.3.5.4.19 Observed analysis: analysis of one full plate of database samples

☐

7.3.5.4.20 Independent analysis: analysis of two batches of previously completed casework known samples

☐

7.3.5.4.21 Independent analysis: analysis of two full plates of previously completed database samples

☐

7.3.5.4.22 Assessment of Analysis section

7.3.5.4.22.1 Competency exam: Analysis of one previously completed casework batch☐7.3.5.4.22.2 Competency exam: Analysis of one previously completed casework knowns batch☐

1)		
2)		

1)		
2)		

--	--

--	--

Approved by Director: Dr. Guy Vallaro

		Trainee Initials/Date	Trainer Initials/Date
7.3.5.4.22.3	<u>Competency exam:</u> Analysis of one previously completed full plate of database samples		
<input type="checkbox"/>			
7.3.5.5	CODIS Training		
7.3.5.5.1	Hands-on training: Deducing CODIS profiles from mixtures using and without using, elimination knowns		
<input type="checkbox"/>			
7.3.5.5.2	SOP Review: DNA SOP-13 CODIS Profile Entry and Data Bank		
<input type="checkbox"/>			
7.3.5.5.3	Independent analysis: Complete worksheet for CODIS eligibility and entry of 5 samples		
<input type="checkbox"/>			
7.3.5.5.4	Independent analysis: Complete required CODIS elements per State CODIS Administrator, including all applicable examinations on CODIS WAN		
<input type="checkbox"/>			
7.3.5.5.5	Written exam: determining eligibility of CODIS samples & profiles to be entered		
<input type="checkbox"/>			
7.3.5.5.6	DNA QR-284, CODIS Training, completed		
<input type="checkbox"/>			
7.3.5.6	Analytical Comparisons, STRmix and Report Writing		
	This section must be completed by Forensic Science Examiners 1, 2, and 3, as well as CCTs to become Forensic Science Examiner 1s.		
7.3.5.6.1	Goals:		
	To provide training in interpretation of DNA results (including mixtures, single source profiles, parentage testing); to write a comprehensive report that accurately reflects the DNA typing results, to learn case flow from beginning to end.		

Approved by Director: Dr. Guy Vallaro

7.3.5.6.2	Tasks	Trainee Initials/Date	Trainer Initials/Date
7.3.5.6.2.1 <input type="checkbox"/>	SOP Review: DNA SOP-6 Report Templates, Statistic Templates, Stutters, and Criticals.		
7.3.5.6.2.2 <input type="checkbox"/>	Review case flow operations through current DNA WIs		
7.3.5.6.2.3 <input type="checkbox"/>	Reading: Forensic DNA Evidence Interpretation Second Edition, Chapters 8 & 9		
7.3.5.6.2.4 <input type="checkbox"/>	Reading: Read STRmix Operation Manual		
7.3.5.6.2.5 <input type="checkbox"/>	SOP Review: SOPs and workflows for STRmix		
7.3.5.6.2.6 <input type="checkbox"/>	Review of all PowerPoint presentations contained in in STRmix training folder, completing all associated quizzes and worksheets that complement the presentations (including likelihood ratio training)		
7.3.5.6.2.7 <input type="checkbox"/>	STRmix software demonstration		
7.3.5.6.2.8 <input type="checkbox"/>	Practice: Complete STRmix analysis of sample profiles, completing 2nd diagnostic worksheets and number of contributor assessment		
7.3.5.6.2.9 <input type="checkbox"/>	Practical test: Complete STRmix practical tests #1 and #2, including report writing	1) 2)	
7.3.5.6.2.10 <input type="checkbox"/>	Read GL-11 Control of Records for lab polices on case file confidentiality, contents of file and discussion of lab results to submitting agencies and attorneys.		

Approved by Director: Dr. Guy Vallaro

		Trainee Initials/Date	Trainer Initials/Date
7.3.5.6.2.11	Use of report templates/macros and standardized wording for report consistency		
<input type="checkbox"/>			
7.3.5.6.2.12	Report Workbook macro training		
<input type="checkbox"/>			
7.3.5.6.2.13	Understand writing of written notification, such as amended letters and reports, and change in CODIS profile memos.		
<input type="checkbox"/>			
7.3.5.6.2.14	Observation of report writing, to include a minimum of 20 reports, 10 at least with known comparisons and 5 at least with STRmix		
<input type="checkbox"/>			
7.3.5.6.2.15	Documentation of reports in the LIMS computer system		
<input type="checkbox"/>			
7.3.5.6.2.16	Training on finalizing reports and submitting finalized reports to submitting agencies.		
<input type="checkbox"/>			
7.3.5.6.2.17	Observed review: Technical Review of 10 cases		
<input type="checkbox"/>			
7.3.5.6.2.18	Observed review: Administrative Review of 3 cases		
<input type="checkbox"/>			
7.3.5.6.2.19	Independent work: Write reports on one mock batch set for a minimum of 5 reports, to be inclusive of all aspects of Nuclear DNA report writing. Documented on QR-4s kept in training binder along with notes for corrections.		
<input type="checkbox"/>			
7.3.5.6.2.20	Competency exam: One previously analyzed mock batch (3 cases minimum) report writing, to include comparisons and STRmix		
<input type="checkbox"/>			

Approved by Director: Dr. Guy Vallaro

		Trainee Initials/Date	Trainer Initials/Date
<input type="checkbox"/>	7.3.5.6.2.21 Moot Court, to include all aspects of analysis		
<input type="checkbox"/>	7.3.5.6.3 Parentage		
<input type="checkbox"/>	7.3.5.6.3.1 Review SOP pertaining to Criminal Parentage Testing		
<input type="checkbox"/>	7.3.5.6.3.2 Hands-on Training: Statistics for Criminal Parentage testing, to include worksheet where independent statistics are manually calculated.		
<input type="checkbox"/>	7.3.5.6.3.3 Hands-on Training: Report writing for criminal parentage testing		
<input type="checkbox"/>	7.3.5.6.3.4 Observed training on parentage report writing and statistics		
<input type="checkbox"/>	7.3.5.6.3.5 Overview of Paternity Calculations, view PowerPoint		
<input type="checkbox"/>	7.3.5.6.3.6 Competency test on parentage comparisons with report writing, written test and oral test (can be moot court, included with others)		
7.3.5.7	Male Screen Analysis/Data Ownership		
<input type="checkbox"/>	7.3.5.7.1 Review of entire male screen process with Trainee		
<input type="checkbox"/>	7.3.5.7.2 SOP overview: SOP-34 and SOP-6		
<input type="checkbox"/>	7.3.5.7.3 Observed analysis: two male screen batches, entire process including creation of LIMS requests (analyst observed by trainee)	1)	
		2)	

Approved by Director: Dr. Guy Vallaro

		Trainee Initials/Date	Trainer Initials/Date
7.3.5.7.4	Independent analysis: one mock male screen batch		
<input type="checkbox"/>			
7.3.5.7.5	Observe analysis: creation of halt and quant and amplification sheets from batch containing male Screen samples		
<input type="checkbox"/>			
7.3.5.7.6	Creation of halt at quant and amp sheet for two previously reported batches	1) <input type="text"/>	<input type="text"/>
<input type="checkbox"/>		2) <input type="text"/>	<input type="text"/>
7.3.5.7.7	Competency exam, including mock male screen Batch, creation of halt at quant and amp sheets from mock casework batch, written exam and moot court questions.	<input type="text"/>	<input type="text"/>
<input type="checkbox"/>			
7.3.5.8	Legal Issues		
	This section must be completed by Forensic Science Examiners 1, 2, and 3, laboratory assistants, and CCTs in the DNA section.		
7.3.5.8.1	Goals		
7.3.5.8.1.1	To give the examiner an overview of the criminal justice system regarding expert witness testimony		
7.3.5.8.1.2	To become familiar with the legal requirements for testimony in Connecticut with the expectations of the CT DESPP Division of Scientific Services		
7.3.5.8.2	Tasks		
7.3.5.8.2.1	Receive guidance on examiner qualifications (voir dire), courtroom appearance (attire and demeanor), court structure, and pertinent rules of a courtroom. Make analyst aware to review notes prior to testimony.	<input type="text"/>	<input type="text"/>
<input type="checkbox"/>			

Approved by Director: Dr. Guy Vallaro

		Trainee Initials/Date	Trainer Initials/Date
7.3.5.8.2.2	Prepare a curriculum vitae able to be presented as part of a discovery process.		
<input type="checkbox"/>			
7.3.5.8.2.3	Review discovery and admissibility rules located on S: Drive Legal Training Folder		
<input type="checkbox"/>			
7.3.5.8.2.4	Review ethical responsibility of an expert witness, included in GL 1.4 and GL 5		
<input type="checkbox"/>			
7.3.5.8.2.5	Read DOJ Code of Professional Responsibility for the Practice of Forensic Science		
<input type="checkbox"/>			
7.3.5.8.2.6	Read: Transcripts and testimony for presentation of DNA test results		
<input type="checkbox"/>			
7.3.5.8.2.7	Review Federal DNA Identification act: https://www.govregs.com/uscode/34/12592		
<input type="checkbox"/>			
7.3.5.8.2.8	Read specific chapters of "Strengthening Forensic Science in the United States" (Chapters 1, 3, 4, 5 (only as it pertains to trainee's job scope), & 7)		
<input type="checkbox"/>			
7.3.5.8.2.9	Read: State admissibility standard (State V. Porter)		
<input type="checkbox"/>			
7.3.5.8.2.10	Read: Federal admissibility standards (Frye, Daubert)		
<input type="checkbox"/>			
7.3.5.8.2.11	Read: Compilation of Connecticut Statutes related to DNA		
<input type="checkbox"/>			
7.3.5.8.2.12	Observe expert testimony		
<input type="checkbox"/>			

Approved by Director: Dr. Guy Vallaro

Trainee Initials/Date	Trainer Initials/Date

- 7.3.5.8.2.13 ☐ Participate in a practice testimony, including direct and cross examinations

- 7.4 **Estimated Schedule for Module Completion:** (Based on time for new employee with no previous experience.)

<u>Section #</u>	<u>Estimated Time to Complete</u>
7.3.1-4	4-8 weeks
7.3.5-7	4 months

Technical Leader of the DNA Section will approve all the training. Currently qualified analysts that have completed this training manual may use DNA QR-283 "Addendum Training Record" to document successful completion of new methodologies and/or platforms. This record will be kept with all other training records.

7.5 Authorizations – Nuclear DNA Casework

The analyst will perform casework and/or reviews and after gaining experience in conducting analysis and/or reporting findings in a variety of evidence type/conclusions in the discipline. Their experience will be evaluated for the ability to conduct technical reviews. Once training has been completed, they may be deemed competent to conduct administrative reviews. Once the analyst has completed the required number of cases (see below), they may be deemed competent to conduct technical and/or administrative reviews.

Trainee Initials/Date	Trainer Initials/Date

- 7.5.1 ☐ For Technical Reviews: New Analyst: complete 15 total cases
- 7.5.2 ☐ For Technical Reviews: Analyst with prior experience: complete 8 total cases

7.6 Kinship Analysis

Upon successful completion of the DNA training program select personnel may be trained in Kinship Analysis. This program outlines formal procedures for the training of personnel in Kinship analysis procedures and defines how an individual examiner will be assessed for competency. The goal of this training program is to develop an examiner capable of performing independent Kinship analysis for non-criminal cases. For examiners that have previous experience in Kinship analysis methods, the training program may be abbreviated as warranted with approval of the Technical Leader.

This training program is designed to supplement successful college coursework in biochemistry, molecular biology, genetics, and population genetics/statistics as required by the FBI QAS. The training program will expose the DNA trainee to the scientific principles underlying Kinship analysis used by the Laboratory.

Assessments will be made at the completion of each module. Each module will be approved by the technical leader once the DNA trainee has been deemed competent at the specified tasks. The DNA trainee will maintain a training folder containing the results of each module as they are being completed. The final paperwork will be retained by the Laboratory. DNA Section personnel must pursue continuing education in the form of workshops, seminars, professional conferences or collegiate coursework as well as further on-the-job training as specified by the FBI QAS standards.

The Laboratory complies with the coursework requirements set forth in the Scientific Working Group on DNA Analysis Methods (SWGDAM) *Guidelines for a Quality Assurance Program for DNA Analysis* and the DNA Advisory Board's standards, *Quality Assurance Standards for Forensic DNA Testing Laboratories*. Examiners must have completed coursework and/or training in Molecular Biology, Genetics, Biochemistry, and Population Genetics/Statistics prior to performing casework.

7.6.1 Goal: Upon completion, the examiner will be familiar with kinship calculations utilized for non-criminal cases.

7.6.2 Tasks

7.6.2.1 SOP Review: DNA SOP-25, Kinship Analysis

☐

7.6.2.2 Review: CT DESPP Internal Kinship Validation and summaries

☐

Trainee Initials/Date	Trainer Initials/Date

Approved by Director: Dr. Guy Vallaro

- | | Trainee
Initials/Date | Trainer
Initials/Date |
|---|--------------------------|--------------------------|
| 7.6.2.3
<input type="checkbox"/> Reading: SWGDAM Guidelines for Missing Persons Casework (current version) | | |
| 7.6.2.4
<input type="checkbox"/> Hands-on training: Practice exercises in Popstats | | |
| 7.6.2.5
<input type="checkbox"/> Hands-on training: Perform kinship calculations by hand | | |
| 7.6.2.6
<input type="checkbox"/> Competency: Written evaluation & moot court detailing trainees understanding of kinship analysis and calculations | | |

*Approved by Director: Dr. Guy Vallaro***7.7 Competencies Memos Received**

Method(s) covered in competency memo

Date Received

_____	_____
_____	_____
_____	_____
_____	_____
_____	_____
_____	_____
_____	_____
_____	_____
_____	_____
_____	_____
_____	_____
_____	_____
_____	_____
_____	_____
_____	_____
_____	_____

7.8 Approval of training plan by Technical Leader

The signature of the Technical Leader below signifies his/her review and approval of the individualized training program & tasks designated for Trainee to complete. Training will not commence until after approval of Technical Leader is obtained.

Technical Leader_____
Date

7.9**Approval of Training and Permission to Start Competency Examinations**

This page can be duplicated as needed to account for different stages in the training process in which trainee would require approval for competency exams at different times. The completion of this form supersedes the requirement for the completion of GL 14.3. Trainee can be deemed competent in individual areas prior to final approval.

Analyst: _____

Training for: _____

To include the following competencies:

☐ 7.3.1.4.1☐ 7.3.2.4.1☐ 7.3.3.4.1☐ 7.3.4.4.1☐ 7.3.5.1.2.36.6☐ 7.3.5.1.3.1☐ 7.3.5.1.3.2☐ 7.3.5.1.3.3☐ 7.3.5.2.2.25☐ 7.3.5.3.13☐ 7.3.5.4.22.1☐ 7.3.5.4.22.2☐ 7.3.5.4.22.3☐ 7.3.5.5.5☐ 7.3.5.6.2.20☐ 7.3.5.6.2.21☐ 7.3.5.6.3.6☐ 7.3.5.7.7

Initialing below states that all training materials have been reviewed, and training in area noted above has been completed to satisfaction so that the above competencies can commence.

Initials

Date

Training binder reviewed by: Technical Leader: _____

Assistant/Deputy Director: _____

Quality Manager: _____

*Approved by Director: Dr. Guy Vallaro***7.10 Final approval of Technical Leader**

The signature of the Technical leader below signifies the successful completion of all tasks assigned to trainee in this manual, assigned to such trainee at his/her commencement of employment at the DESPP Division of Scientific Services DNA Unit:

Technical Leader_____
Date**7.11 Assistant Director or Quality Manager Review**_____
Signature_____
Date**7.12 Deputy Director Review**_____
Signature_____
Date