

Acquisition Method Report



Acquisition Method Info

|                    |                                                 |                        |                           |
|--------------------|-------------------------------------------------|------------------------|---------------------------|
| Operator Name      | SYSTEM (SYSTEM)                                 |                        |                           |
| Created By         | SYSTEM (SYSTEM)                                 | Creation Date          | 2024-08-01 10:17:28-04:00 |
| Last Modified By   | SYSTEM (SYSTEM)                                 | Last Modification Date | 2024-08-02 11:56:43-04:00 |
| Method Name        | GHB.m                                           |                        |                           |
| Method Path        | D:\Projects\TOX Lab\Methods\GHB.m (Version: NA) |                        |                           |
| Method Description |                                                 |                        |                           |
| Pre Run Script     |                                                 |                        |                           |
| Post Run Script    |                                                 |                        |                           |
| Est. Run Time      | 5.00                                            |                        |                           |
| Device List        |                                                 |                        |                           |
|                    | Multisampler                                    |                        |                           |
|                    | Binary Pump                                     |                        |                           |
|                    | Column Comp.                                    |                        |                           |
|                    | TQ                                              |                        |                           |

Name: MultisamplerModule: G7167A

Sampling Speed

|                         |              |
|-------------------------|--------------|
| Draw Speed              | 100.0 µL/min |
| Eject Speed             | 400.0 µL/min |
| Wait Time After Drawing | 1.2 s        |

Injection

|                  |               |
|------------------|---------------|
| Needle Wash Mode | Standard Wash |
| Injection Volume | 5.00 µL       |

Standard Needle Wash

|                  |            |
|------------------|------------|
| Needle Wash Mode | Flush Port |
| Duration         | 3 s        |

High Throughput

|                                                      |     |
|------------------------------------------------------|-----|
| Injection Valve to Bypass for Delay Volume Reduction | No  |
| Sample Flush-Out Factor                              | 5.0 |

Overlapped Injection

|                           |    |
|---------------------------|----|
| Overlap Injection Enabled | No |
|---------------------------|----|

Needle Height Position

|                              |        |
|------------------------------|--------|
| Draw Position Offset         | 0.0 mm |
| Use Vial/Well Bottom Sensing | No     |

Stoptime

|               |          |
|---------------|----------|
| Stoptime Mode | Time set |
| Stoptime      | 4.00 min |

Posttime

# Acquisition Method Report



Posttime Mode Off

## Use Injector Program

Use Injector Program No

Name: Binary Pump

Module: G7112B

Flow 0.300 mL/min  
Use Solvent Types Yes  
Low Pressure Limit 0.00 bar  
High Pressure Limit 600.00 bar  
Maximum Flow Gradient 100.000 mL/min<sup>2</sup>

## Stroke A

Automatic Stroke Calculation A Yes

## Stroke B

Automatic Stroke Calculation B Yes

## Stoptime

Stoptime Mode Time set

Stoptime 4.00 min

## Posttime

Posttime Mode Off

## Solvent Composition

|   | Channel | Solvent 1 | Name 1 | Solvent 2 | Name 2 | Selected | Used | Percent (%) |
|---|---------|-----------|--------|-----------|--------|----------|------|-------------|
| 1 | A       | H2O       |        | H2O       |        | Ch. 2    | Yes  | 95.0 %      |
| 2 | B       | H2O       |        | MeOH      |        | Ch. 2    | Yes  | 5.0 %       |

## Timetable

|   | Time (min)       | A (%)  | B (%) | Flow (mL/min) | Pressure (bar) |
|---|------------------|--------|-------|---------------|----------------|
| 1 | Start. Cond. min | 95.0 % | 5.0 % | 0.300 mL/min  | 600.00 bar     |
| 2 | 3.00 min         | 92.0 % | 8.0 % | 0.300 mL/min  | 600.00 bar     |
| 3 | 3.01 min         | 95.0 % | 5.0 % | 0.300 mL/min  | 600.00 bar     |
| 4 | 4.00 min         | 95.0 % | 5.0 % | 0.300 mL/min  | 600.00 bar     |

Name: Column Comp.

Module: G7116A

## Left Temperature Control

Left Temp. Control Mode Temperature Set

Left Temperature 30.0 °C

## Enable Analysis Left Temperature

Enable Analysis Left Temperature On Yes

Enable Analysis Left Temperature Value 0.8 °C

Left Temp. Equilibration Time 0.0 min

## Right Temperature Control

Right temperature Control Mode Combined

# Acquisition Method Report



## Enable Analysis Right Temperature

Enable Analysis Right Temperature On Yes  
 Enable Analysis Right Temperature Value 0.8 °C  
 Right Temp. Equilibration Time 0.0 min

## Enforce column for run

Enforce column for run is enabled No

## Thermal equilibration devices

Enforce left thermal equilibration device No  
 Enforce right thermal equilibration device No

## Stoptime

Stoptime Mode As Pump/Injector

## Posttime

Posttime Mode Off

## Timetable

Valve Position Position 2 (Port 1 -> 2)  
 Position Switch After Run Do not switch

Name: TQ Mass Spectrometer

Module: G6475A

## General

Tune File atunes.tune  
 Ion Source AJS ESI  
 Optimized Method Off  
 Last Compound Optimized Time  
 Last Source Optimized Time  
 Time Filter Enabled Off

## Stop Time

Stop Mode By Pump Time

## Time Segment

|   | Start Time (min) | Scan Type | Detector Gain Factor (+) | Detector Gain Factor (-) | Store Data |
|---|------------------|-----------|--------------------------|--------------------------|------------|
| 1 | 0                | MRM       | 1                        | 1                        | Yes        |

## Time Segment 1

Scan Type: MRM

|   | Compound group | Compound name | Compound formula | Ion species | CAS | z | Monoisotopic mass (m/z) | ISTD? | Precursor (m/z) | MS1 res | Product (m/z) | MS2 res | Dwell (ms) | Fragment (V) | CAV (V) | CE (V) | Polarity |
|---|----------------|---------------|------------------|-------------|-----|---|-------------------------|-------|-----------------|---------|---------------|---------|------------|--------------|---------|--------|----------|
| 1 |                | GHB-D6        | C4HD6O3          | (M+H)+      |     | 1 | 109.1 m/z               | No    | 110.1 m/z       | Unit    | 92 m/z        | Unit    | 100 ms     | 60 V         | 4 V     | 5 V    | +        |
| 2 |                | GHB-D6        | C4HD6O3          | (M+H)+      |     | 1 | 109.1 m/z               | No    | 110.1 m/z       | Unit    | 49 m/z        | Unit    | 100 ms     | 60 V         | 4 V     | 20 V   | +        |
| 3 |                | GHB           | C4H8O3           | (M+H)+      |     | 1 | 104 m/z                 | No    | 105 m/z         | Unit    | 87 m/z        | Unit    | 100 ms     | 60 V         | 4 V     | 5 V    | +        |
| 4 |                | GHB           | C4H8O3           | (M+H)+      |     | 1 | 104 m/z                 | No    | 105 m/z         | Unit    | 45 m/z        | Unit    | 100 ms     | 60 V         | 4 V     | 20 V   | +        |

# Acquisition Method Report



|   | Compound group | Compound name | Compound formula | Ion species | CAS | z | Monoisotopic mass (m/z) | ISTD? | Precursor (m/z) | MS1 res | Product (m/z) | MS2 res | Dwell (ms) | Fragmentor (V) | CAV (V) | CE (V) | Polarity |
|---|----------------|---------------|------------------|-------------|-----|---|-------------------------|-------|-----------------|---------|---------------|---------|------------|----------------|---------|--------|----------|
| 5 |                | GBL           | C4H6O2           | (M+H)+      |     | 1 | 86 m/z                  | No    | 87 m/z          | Unit    | 45 m/z        | Unit    | 100 ms     | 70 V           | 4 V     | 15 V   | +        |
| 6 |                | GBL           | C4H6O2           | (M+H)+      |     | 1 | 86 m/z                  | No    | 87 m/z          | Unit    | 43 m/z        | Unit    | 100 ms     | 70 V           | 4 V     | 10 V   | +        |

## Source Parameter

|   | Parameter Id                | Positive Value | Negative Value |
|---|-----------------------------|----------------|----------------|
| 1 | Gas Temperature (°C)        | 300            | 300            |
| 2 | Gas Flow (L/min)            | 13             | 13             |
| 3 | Nebulizer (psi)             | 35             | 35             |
| 4 | Capillary Voltage (V)       | 4000           | 4000           |
| 5 | Sheath Gas Temperature (°C) | 250            | 250            |
| 6 | Sheath Gas Flow (L/min)     | 11             | 11             |
| 7 | Nozzle Voltage (V)          | 1500           | 1500           |

## Optimizer Settings

### Optimizer Compound Settings

|                                 |          |
|---------------------------------|----------|
| Precursor abundance threshold   | 5000     |
| Product ion abundance threshold | 500      |
| MRM abundance threshold         | 100      |
| Optimizer Dwell Time            | 130 ms   |
| Maximum number of product ions  | 4        |
| Product ion low m/z cut-off     | Disabled |

|   | Name             | Option                               | From | To  | Step Size |
|---|------------------|--------------------------------------|------|-----|-----------|
| 1 | Fragmentor (V)   | Find product ions from defined range | 110  | 140 | 10        |
| 2 | Collision Energy | Find product ions from defined range | 0    | 20  | 4         |

## Chromatogram

|   | Chromatogram Type | Label |
|---|-------------------|-------|
| 1 | TIC               | TIC   |

## Timetable

# Acquisition Method Report



|   | Start Time (min) | Timetable Type | Timetable Value |
|---|------------------|----------------|-----------------|
| 1 | 0                | Diverter Valve | To Waste        |
| 2 | 1.4              | Diverter Valve | To MS           |
| 3 | 2.6              | Diverter Valve | To Waste        |

## Advanced override settings

### Positive

|   | Parameter | Override Tune | Dynamic Mass (m/z) | Tune Setting | Override Value |
|---|-----------|---------------|--------------------|--------------|----------------|
| 1 | HED (kV)  | No            | m/z                | Standard     | Large molecule |

### Negative

Name: DA - Intelligent Reflex Parameters

Intelligent Reflex Enabled

No