

Document Title: Screening Tests for Blood

Controlled: Yes, with red stamp present

Controlled By: Quality Manager

Prepared By: \_\_\_\_\_ Date: \_\_\_\_\_

Approved By: \_\_\_\_\_ Date: \_\_\_\_\_

**A. PURPOSE:**

To perform screening tests for the presence of blood in Forensic samples.

**B. RESPONSIBILITY:**

Forensic Science Examiners from the Connecticut State Forensic Science Laboratory who have been trained in the discipline of blood screening according to SOP-FB-31 (Training Manual).

**C. SAFETY:**

Use appropriate measures for the proper handling of the o-Tolidine solution according to SOP-GL-2 (Safety Manual).

**D. DEFINITIONS:**

1. PBS: Phosphate Buffered Saline
2. KM: Kastle-Meyer Test

**E. PROCEDURE:**

These tests will be performed at the discretion of the examiner based on the submitting agency requests, case information and the condition of the evidence.

1. Materials:
  - a. Phenolphthalin solution (KM test)
  - b. o-Tolidine solution
  - c. 3.0% Hydrogen peroxide
  - d. Controls: positive (known blood stain) and negative (blank filter paper)
  - e. PBS
  - f. 0.5% ammonia
  - g. Cotton swabs or spot plates
2. Procedure:
  - a. Test a positive control and a negative control (blank filter paper) with PBS according to the following procedure (steps 2.b. – 2.e.) prior to the questioned sample.
    - aa. If controls yield the appropriate results, record on the appropriate Quality Record Worksheet and test the questioned samples.

- E. 2. a. bb. If controls do not yield the appropriate results, review the procedure and retest the controls prior to the questioned sample.
  - b. Moisten swab with PBS to remove and test a portion of the questioned sample or stain. A portion of the questioned sample or stain may be removed and tested directly.
    - aa. 0.5% ammonia may be used in place of PBS for aged samples.
    - bb. Both controls must be tested with 0.5% ammonia according to steps 2.a. – 2.e., prior to testing the questioned sample.
    - cc. If 0.5% ammonia is used, record on the appropriate Quality Record Worksheet and the General Reagent Sheet (FBQR-09).
  - c. Add one drop of phenolphthalin or o-Tolidine solution to the portion of questioned sample. *Caution:* Testing with o-Tolidine solution should be performed under the hood.
  - d. If no color change occurs, add one (1) drop of 3% H<sub>2</sub>O<sub>2</sub>.
  - e. Observe any color change within 10 seconds.
3. Results:
- a. *Positive.* If blood or "other peroxidase-type material" is present, a color change will occur after the addition of 3% H<sub>2</sub>O<sub>2</sub>. Phenolphthalin (KM test) will turn pink; o-Tolidine will turn blue.
  - b. *Negative.* A negative reaction will show no color change within 10 seconds of the addition of 3% H<sub>2</sub>O<sub>2</sub>, indicating no blood was detected.
  - c. *Inconclusive.* The appearance of a color change without the addition of 3% H<sub>2</sub>O<sub>2</sub> may indicate the presence of a chemical oxidant in the stain.
  - d. Record the results on the appropriate Quality Record Worksheet.
  - e. If the quantity of stain is insufficient to perform a confirmatory test, a positive chemical screening test is sufficient to forward the stain directly for DNA analysis.
4. Record reagent(s) used on the General Reagent Sheet (FBQR-09).

**F. REFERENCES:**

1. Holland, V.R.B., C.Saunders, F.L.Rose and A.L.Wulpoll. A safer substitute for benzidine in the detection of blood. Tetrahedron, Vol.30, 1974, pp. 3299-3302.
2. Hunt, A.C., C. Corley, B.E. Dodd and F.E. Camps. The identification of human blood stains: critical survey. J.Forensic Med. Vol. 7, 1960, pp. 112-130.
3. Koan, J. and T. O'Kelly. An ortho-tolidine method for the detection of occult blood in feces. J. Clinical Pathol. Vol.8, 1955, pp. 249-251.
4. Ruttan,R.F. and R.H.M. Hardisty. A new reagent for detecting occult blood. Can. Med. Associat. J. Vol.41(n.s.2), 1912, pp. 995-998.
5. Metropolitan Police Forensic Science Laboratory. Biology Methods Manual. 1978, pp. 2-88 to 2-90.
6. SOP-GL-2 (Safety Manual).