

Document Title: Species Double Diffusion Test (Ouchterlony)

Controlled: Yes, with red stamp present

Controlled By: Quality Manager

Prepared By: _____ Date: _____

Approved By: _____ Date: _____

A. PURPOSE:

To determine the species of origin in cases when animal blood is suspected.

B. RESPONSIBILITY:

Forensic Science Examiners from the Connecticut State Forensic Science Laboratory who have been trained in the discipline of the species double diffusion test procedure according to SOP-FB-31 (Training Manual).

C. SAFETY:

Use appropriate measures for the proper handling of the Ouchterlony plates according to SOP-GL-2 (Safety Manual).

D. DEFINITIONS:

PBS: Phosphate Buffered Saline

E. PROCEDURE:

This procedure will be performed at the discretion of the examiner based on the submitting agency requests, case information and the condition of the evidence.

1. Materials:

- a. Ouchterlony plates
- b. Animal anti-sera (anti-dog, -cat, -deer, etc.)
- c. Known bloodstain extracts or thawed sera controls
- d. PBS
- e. 0.5% ammonia (NH₄)
- f. Disposable pipets or micropipet and tips
- g. Microcentrifuge tubes

2. Procedure:

a. Preparing extracts:

- aa. Extract a portion of the sample in a microcentrifuge tube with 30-50µl of PBS. If necessary, sample can be extracted for a longer period of time on a shaker at room temperature or overnight at 4°C.

- E. 2. a. bb. 0.5% ammonia may be used in place of PBS to extract aged samples or samples that are difficult to extract.
- b. Each plate must contain a positive control and negative reagent blank control.
- aa. Known blood stain extracts, used as positive controls for animal anti-sera, should be made up fresh for each use with the same extraction solution used for the questioned samples. Store remaining extract(s) in the refrigerator and discard after six (6) months or sooner if necessary.
- bb. Frozen aliquots of animal sera may be thawed and used as positive controls for animal anti-sera. Store thawed, remaining animal sera in the refrigerator and discard after six (6) months or sooner if necessary.
- cc. PBS or 0.5% NH₄ may be used as a negative reagent blank control according to the extraction solution used.
- dd. Record the extraction solution(s) used and the animal standards used, if applicable, on the Ouchterlony Quality Record Worksheet (FBQR-08). Record the extraction solutions used for questioned samples on the General Reagent Sheet (FBQR-09).
- c. Show placement of the controls, samples and anti-serum/anti-sera on the schematic diagram of the Ouchterlony Quality Record Worksheet (FBQR-08).
- d. Fill the wells in the Ouchterlony plate with a pipet according to the diagram. Avoid air bubbles or over flowing the wells.
- e. Allow plate to sit upright at room temperature until all samples have diffused into the plate.
- f. Turn the plate upside down and allow the plate to stand overnight (12-16 hrs) at 4°C.
- g. Examine the plate for precipitin lines using back lighting and record the results on the schematic diagram of the Ouchterlony Quality Record Worksheet (FBQR-08).
- Include a copy of this worksheet in the case jacket, file the original in the designated notebook.
- h. If no precipitin line is observed, re-examine the plate for up to 48 hours.
- i. If the controls do not give the appropriate results, the test is considered inconclusive and the sample(s) should be re-run.

E. 3. Results:

- a. *Positive.* A white precipitin line between a sample well and the anti-serum well indicates a positive result.
- b. *Negative.* No line will appear between a sample well and the anti-serum well, indicating a negative result.
- c. *Inconclusive.*
 - aa. No distinguishable white precipitin line observed and/or insufficient extraction of sample.
 - bb. The controls do not give the appropriate results and there is not enough sample to repeat the test.
- d. Record the results on the appropriate Quality Record Worksheet. Note: The reason a result is determined to be inconclusive must also be recorded.

F. REFERENCES:

1. Ouchterlony, O., 1948a, "Antigen-antibody reactions in gels," Acta. Pathol. Microbiol. Scand. 26 (1949), 507.
2. Ouchterlony, O. "Antigen-antibody reaction in gels", Ark. Kemi. Mineral Geol. 26B (14).
3. Ouchterlony, O. 1949b. Antigen-antibody reactions in gels II. Factors determining the site of the precipitate. Ark. Kemi. 1:43-48.
4. Ouchterlony, O. 1949c. Antigen-antibody reactions in gels III. The time factor. Ark. Kemi. 1:55-59.
5. Ouchterlony, O. 1968. Handbook of Immunodiffusion and Immuno-electrophoresis, Ann Arbor Science Publishers, Inc., Ann Arbor, Michigan.
6. SOP-GL-2 (Safety Manual).