

**Connecticut Department of Public Safety  
Division of Scientific Services  
Forensic Laboratory**

Document ID: SOP-FB-14  
Revision #: 0  
Revision Date: 01/01/2011

Page 1 of 4

Document Title: Rapid Immunoassay Tests for Human Semen  
Controlled: Yes, with red stamp present  
Controlled By: Quality Manager

Prepared By: \_\_\_\_\_ Date: \_\_\_\_\_

Approved By: \_\_\_\_\_ Date: \_\_\_\_\_

**A. PURPOSE:**

To identify Semenogelin (RSID™-Semen) or p30 (ABAcad®) in Forensic samples.

**B. RESPONSIBILITY:**

Forensic Science Examiners from the Connecticut State Forensic Science Laboratory who have been trained in the discipline of rapid immunoassay tests for human semen according to SOP-FB-31 (Training Manual).

**C. DEFINITIONS:**

1. RSID™: Rapid Stain Identification
2. ABAcad®: Rapid Immunoassay

**D. PROCEDURE:**

RSID™-Semen and/or p30 ABAcad® tests will be performed at the discretion of the examiner based on the submitting agency requests, case information and the condition of the evidence.

1. Materials:
  - a. Extracted samples
  - b. RSID™-Semen test cassettes
  - c. RSID™-Semen running buffer
  - d. ABAcad® test device and enclosed dropper
  - e. Microcentrifuge tubes
  - f. Spot plates
  - g. Micropipet and tips

2. Procedures:

RSID™-Semen Test with Universal Buffer

- a. Bring extracted samples to room temperature before use (approximately 10 minutes).
- b. Label RSID™-Semen cassettes with case and item number.
- d. Using a micropipet, add 100µl of extract to well 'S' of the cassette. Note the time immediately after adding the sample.

- D. 2. e. Monitor progress of test results for a 10 minute period. Record final result at 10 minutes. DO NOT record any changes which occur after 10 minutes. Any change in test results after 10 minutes is invalid.
- f. Subsequently, the ABACard® p30 procedure can be run with the remaining extract.

RSID™-Semen Test with Extraction Buffer or PBS


- a. Bring extracted samples to room temperature before use (approximately 10 minutes).
- b. Label RSID™-Semen cassettes with case and item number.
- c. Using a micropipet, place 20µl of the supernatant into a microcentrifuge tube or spot plate and add 80µl of RSID™-Semen running buffer for a total extract volume of 100µl.
- d. Add 100µl of extract to well 'S' of the cassette. Note the time immediately after adding the sample.
- e. Monitor progress of test results for a 10 minute period. Record final result at 10 minutes. DO NOT record any changes which occur after 10 minutes. Any change in test results after 10 minutes is invalid.
- f. Subsequently, the ABACard® p30 procedure can be run with the remaining extract.

ABACard® p30 Test

- a. Bring extracted samples and test devices to room temperature before use (approximately 10 minutes).
- b. Label ABACard® with case and item numbers.
- c. Using the enclosed dropper, add 5-6 drops (~150-200µl) of supernatant to well "S" of the test device. Note the time immediately after adding the sample.
- d. Monitor progress of test results for a 10-minute period. Record final result at 10 minutes. DO NOT record any changes which occur after 10 minutes. Any change in test results after 10 minutes is invalid.


D 3. Results and Interpretation:

RSID™-Semen Test



- a. *Negative.* A visible red line at the Control 'C' position only, indicates a negative result.  
*No human semen detected or the sample is below the detectable level.*
- aa. "High Dose Hook Effect" refers to weak a positive or false negative result seen when very high levels of human semen are present in the sample.
- bb. In such cases or when a strong Acid Phosphatase reaction is observed, the extract may be diluted 1:20 and rerun with the RSID™-Semen Extraction Buffer or Universal Buffer procedure.
- b. *Positive.* Visible red lines at both the Control 'C' and Test 'T' positions indicate a positive result. *Human semen detected.*
- c. *Failed(Invalid).* No visible red line at the Control 'C' position indicates a failed test.  
*No conclusion possible.* Review the test procedure carefully and repeat the test with a new plate.
- d. If a negative result is obtained, the ABACard® p30 procedure below may be used.

ABACard® p30 Test



D. 3. ABACard® p30 Test

- a. *Positive.* If there are two pink lines, one each in the test area 'T' and in the

control area 'C', the test result is *positive* and indicates that the p30 level is present at or above 4ng/ml.

- b. *Negative*. If there is only one pink line in the control area 'C', the test result is *negative*. This may indicate that no p30 is present, the p30 level is below 4ng/ml or the presence of the "High Dose Hook Effect".
  - aa. "High Dose Hook Effect" refers to weak positive or false negative results due to the presence of a high concentration of p30 in the sample.
  - bb. In such cases or when a strong Acid Phosphatase reaction is observed, the extract may be diluted 1:20 and rerun with the ABACard® Procedure.
- c. *Invalid*. If there is no pink line visible in the control area 'C', the test is *inconclusive*. Review the test procedure carefully and repeat the test with a new plate.

4. Record the results of the control(s) and sample(s) on the appropriate Quality Record Worksheet.

#### **E. REFERENCES:**

1. Abacus Diagnostics' *OneStep* ABACard® p30 Test For The Forensic Identification of Semen, provided Technical Information and Protocol sheet.
2. Connecticut State Forensic Science Laboratory, ABACard® p30 Internal Validation, 1998.
3. Independent Forensics, Rapid Stain Identification of Human Semen (RSID™-Semen) Technical Information and Protocol sheets.
4. Old, Dr. Jennifer, Schweers, Dr. Brett A., Boonlayangoor, Dr. P. W., Reich, Dr. Karl, Developmental Validation Studies of RSID™-Semen Lateral Flow Immunochromatographic Strip test for the forensic detection of Seminal Fluid, p 1-36.
5. Connecticut State Forensic Science Laboratory, RSID™ - Semen Internal Validation, 2010.
6. Independent Forensics, RSID™ - Universal Buffer Technical Information and Protocol sheet.
7. Connecticut State Forensic Science Laboratory, RSID™ - Universal Buffer Internal Validation, 2011.