

*Approved by Director: Dr. Guy Vallaro***A. PURPOSE:**

To identify Tamm-Horsfall Protein (RSID™-Urine) in Forensic samples, which indicates the presence of urine.

**B. RESPONSIBILITY:**

Forensic Science Examiners from the Connecticut State Forensic Science Laboratory who have been trained in the discipline of the rapid immunoassay test for urine according to SOP-FB-31 (Training Manual).

**C. DEFINITIONS:**

1. RSID™: Rapid Stain Identification
2. THP: Tamm-Horsfall Protein

**D. PROCEDURE:**

This test will be performed at the discretion of the examiner based on the submitting agency requests, case information and the condition of the evidence.

1. Materials:
  - a. RSID™-Urine Buffer
  - b. RSID™-Urine test cassettes
  - c. Microcentrifuge tubes and spin baskets
  - d. Micropipet and tips
  - e. Wooden sticks
  - f. Ultrasonic bath
  - g. Shaker
  - h. Centrifuge
2. Procedure:
  - a. Do not use this test with samples containing fecal material, refer to SOP-FB-17 (Test for Urobilinogen).


Undiluted urine should not be used with this test, a stain of the sample should be made instead and dried overnight.

- b. Extract a portion of the questioned sample or stain in a microcentrifuge tube with 150µl of RSID™-Urine Buffer and mix with a wooden stick.

*Approved by Director: Dr. Guy Vallaro*

- D. 2. c. Place the sample in an ultrasonic bath for fifteen (15) minutes and then place on the shaker for one (1) hour and forty-five (45) minutes at room temperature (total extraction time of two (2) hours).
- d. If necessary, the sample may be extracted overnight at 4°C after the sonication step.
- e. If the sample was refrigerated, bring to room temperature for approximately ten (10) minutes before testing.
- f. Place the sample in a spin basket and centrifuge for approximately ten (10) minutes within the range of 10,000 – 14,000rpm (13,000 rpm is recommended).
- g. Label the RSID™ cassettes with case and item numbers.
- h. Using a micropipet, add 100µl of extract to well ‘S’ of the cassette. Note the time immediately after adding the sample.
- i. Monitor progress of test results for a fifteen (15) minute period. Record final result at fifteen (15) minutes. **DO NOT** record any changes that occur after fifteen (15) minutes. Any change in the test results after fifteen (15) minutes is invalid.

### 3. Results:



- a. *Negative.* A visible blue line at the Control ‘C’ position only, indicates a *negative* result.  
*No Tamm-Horsfall Protein detected.*
- aa. “High Dose Hook Effect” refers to weak positive or false negative results due to the presence of a high concentration of THP in the sample.
- bb. Under standard laboratory testing and relevant urine concentration ranges, the "High

*Approved by Director: Dr. Guy Vallaro*

Dose Hook Effect" is not observed with the RSID™-Urine Test.

- b. *Positive.* Visible blue lines at both the Control 'C' and Test 'T' positions indicate a positive result. *Tamm-Horsfall Protein detected.*

D. 3. c. *Inconclusive/Failed (Invalid).*

- aa. No visible blue line at the Control 'C' position indicates a failed test. *No conclusion possible.* Review the test procedure carefully and repeat the test with a new plate.

- bb. No distinguishable blue line at the Test 'T' position.

4. Record the results on the appropriate Quality Record Worksheet. Note: The reason a result is determined to be inconclusive must also be recorded.

5. Record test(s) used on the General Reagent Sheet (FBQR-09).

**E. REFERENCES:**

1. Independent Forensics, Rapid Stain Identification of Urine (RSID™ - Urine) Technical Information and Protocol sheet.
2. Old, Dr. Jennifer, Reich, Dr. Karl, Developmental Validation of RSID™ - Urine, p1-19.
3. Connecticut State Forensic Science Laboratory, RSID™ - Urine Internal Validation, 2012.