

Document Title: Test for Creatinine (Jaffe)

Controlled: Yes, with red stamp present

Controlled By: Quality Manager

Prepared By: \_\_\_\_\_ Date: \_\_\_\_\_

Approved By: \_\_\_\_\_ Date: \_\_\_\_\_

**A. PURPOSE:**

To determine the presence of creatinine in a Forensic sample, which indicates the presence of urine.

**B. RESPONSIBILITY:**

Forensic Science Examiners from the Connecticut State Forensic Science Laboratory who have been trained in the discipline of testing for creatinine according to SOP-FB-31 (Training Manual).

**C. SAFETY:**

Use appropriate measures for the proper handling of picric acid solution, sodium hydroxide and glacial acetic acid according to SOP-GL-2 (Safety Manual).

**D. PROCEDURE:**

This test will be performed at the discretion of the examiner based on the submitting agency requests, case information and the condition of the evidence.

1. Materials:

- a. Picric Acid Solution (saturated)
- b. 10% NaOH
- c. Distilled water (dH<sub>2</sub>O)
- c. Glacial Acetic Acid
- d. Controls: positive (known urine stain) and negative (blank filter paper), include substrate control as needed
- e. Disposable pipet or micropipet and tips
- f. Test tubes

2. Procedure:

- a. Test a positive and negative control with the following procedure (steps 2.b. – 2.g.).
  - aa. The controls may be run concurrently with the questioned samples.
  - bb. If limited questioned sample is available, run the controls prior to testing the questioned sample. If controls yield the appropriate results then test the questioned sample.
  - cc. If controls do not yield the appropriate results, review the procedure and retest the controls prior to the questioned samples.

- D. 2. b. Extract a portion of the questioned sample or stain in a test tube or spot plate with enough dH<sub>2</sub>O to cover the sample. Let stand for a minimum of 15 minutes or longer as necessary.
- c. Remove substrate from test tube.
- d. Add equal amounts of picric acid solution and 10% NaOH to an equal amount of the extract.
- e. Allow to stand for 15 minutes and observe the color of the extract.
- f. Add glacial acetic acid in excess (double the total volume of extract, picric acid and 10% NaOH) and observe the color of the extract.
- g. It should be noted that the detection of a urine or ammonia-type odor may also be used in the interpretation of the results (sample may need to be heated at 37°).

3. Results:

- a. *Positive.* A positive result is indicated by a color change from yellow to orange/reddish-orange and develops fully within 15 minutes. A color change back to yellow after the addition of glacial acetic acid, confirms the presence of creatinine.
- b. *Negative.* A negative result is indicated by the bright yellow color of the picric acid.
- c. *Inconclusive.* A sample that remains unchanged after the addition of glacial acetic acid is considered inconclusive.
- d. It is important to compare results against the positive and negative controls.
- e. Record the results of the controls and samples on the appropriate Quality Record Worksheet.
- f. A 2<sup>nd</sup> examiner will observe and confirm results and initial the appropriate Quality Record Worksheet.

4. Record solutions used on the General Reagent Sheet (FBQR-09).

**E. REFERENCES:**

- 1. Jaffe, M. "Uebeden Niederschlag, Welchen Pikrinsäure in normalen Harn erzeugt und über line neue Reaktion des Kreatinins." Z. Physical Chem, 10: 391-400. 1886.
- 2. Metropolitan Police Forensic Science Laboratory. Biology Methods Manual. 1978, pp. 4-4 to 4-5.
- 3. SOP-GL-2 (Safety Manual).