

*Approved by Director: Dr. Guy Vallaro***A. PURPOSE:**

To prepare reagents for blood screening tests and to perform quality control on prepared reagents.

B. RESPONSIBILITY:

Forensic Science Examiners 1 and 2 in the Forensic Biology Section. Ordering information is maintained in a log book in the Forensic Biology Section.

C. SAFETY:

Use appropriate measures for the proper handling of o-Tolidine, ethanol, glacial acetic acid and potassium hydroxide according to SOP-GL-2 (Safety Manual) and the Material Safety Data Sheets.

D. DEFINITIONS:

1. dH₂O: Distilled water
2. KM: Kastle-Meyer Test

E. PROCEDURE:**o-Tolidine Solution****1. Materials:**

- | | | |
|----|---------------------|--------|
| a. | o-Tolidine | 1.6g |
| b. | Ethanol | 40.0ml |
| c. | Glacial acetic acid | 30.0ml |
| d. | dH ₂ O | 30.0ml |

2. Procedure:

- a. Mix all materials together and place into a brown dropper bottle.
- b. Test each new batch of reagent before use according to SOP-FB-07 (Screening Tests for Blood) and the o-Tolidine Reagent Log Sheet. Record the required information.
- c. If the appropriate results are not obtained, discard the reagent, review the procedure and make new reagent.
- d. If the reagent is suitable for use, record the solution, lot # (date of preparation), control date and examiner's initials on the dropper bottle and store in the refrigerator.

Approved by Director: Dr. Guy Vallaro

- E. 2. e. Discard after one (1) year or sooner if the reagent changes from transparent beige to a dark brown color.

Phenolphthalin Stock Solution

1. Materials:

- | | | |
|----|---------------------|-------|
| a. | Phenolphthalin | 2g |
| b. | Potassium hydroxide | 20g |
| c. | dH ₂ O | 100ml |
| d. | Granular Zinc | 20g |

2. Procedure:

- Dissolve potassium hydroxide in dH₂O. Note: The solution will be warm from the reaction.
- Add granular zinc followed by phenolphthalin. Solution turns pink immediately.
- Swirl solution until it becomes colorless.
- Place stock solution in a brown bottle containing zinc.
- Discard after six (6) months or sooner if the reagent turns from colorless to pink.

Phenolphthalin Working Solution

1. Materials:

- | | | |
|----|-------------------------------|---------|
| a. | Phenolphthalin stock solution | 1 part |
| b. | Ethanol | 4 parts |
| c. | Granular Zinc | ~15g |

2. Procedure:

- Dilute stock solution 1:5 in ethanol (1 part to 4 parts).
- Place working solution in brown dropper bottles containing zinc.
- Test each new batch of the working solution before use according to SOP-FB-07 (Screening Tests for Blood) and the KM Reagent Log Sheet. Record the required information.
- If the appropriate results are not obtained, discard the reagent, review the procedure and make new reagent.

Approved by Director: Dr. Guy Vallaro

- e. If the reagent is suitable for use, record the solution, lot # (date of preparation), control date and examiner's initials on the stock and dropper bottles and store in the refrigerator.
- f. Discard after six (6) months or sooner if the reagent turns from colorless to pink.

E. 3 % Hydrogen Peroxide

- 1. This chemical is purchased from an outside vendor and tested/used as received.
- 2. Test the new manufacturer's lot before use according to SOP-FB-07 (Screening Tests for Blood) and the KM and o-Tolidine Reagent Log Sheets. Record the required information.
- 3. If the appropriate results are not obtained, review the procedure, repeat the test and replace the chemical if necessary.
- 4. If the lot is suitable for use, record the date received, date opened and examiner's initials on the stock bottles and store in the refrigerator.
- 5. Record the chemical, lot #, manufacturer's expiration date, fill date and initials on the dropper bottles.
- 6. Discard according to the manufacturer's expiration date or sooner if a decrease in reaction activity is noted.

F. REFERENCES:

- 1. Metropolitan Police Forensic Science Laboratory. Biology Methods Manual. 1978, pp. 2-88 to 2-90.
- 2. SOP-GL-2 (Safety Manual).
- 3. Material Safety Data Sheets.