

**BODY FLUID STANDARDS****22.1 PURPOSE**

To prepare body fluid standards for the purpose of maintaining quality control of reagents and alternate light sources.

**22.2 RESPONSIBILITY**

Forensic Science Examiners in the Forensic Biology Unit.

**22.3 SAFETY**

Use appropriate measures for the proper handling of biohazardous materials according to the GL-2 (Safety Manual).

**22.4 DEFINITIONS**

- A. PTT: Purple Top Tube
- B. KM: Kastle Meyer Test
- C. AP: Acid Phosphatase
- D. ALS: Alternate Light Source
- E. RSID™: Rapid Stain Identification
- F. ABACard®: Rapid Immunoassay
- G. PBS: Phosphate Buffered Saline

**22.5 PROCEDURE**

For all standards prepared, record the appropriate information on the Body Fluid Standard Reagent Log Sheet.

**21.5.1: Materials**

- A. Body fluid samples (human unless otherwise specified)
- B. PBS
- C. dH<sub>2</sub>O
- D. Cloth swatches (white and black)
- E. Filter paper
- F. Swabs
- G. Glass slides
- H. Coin envelopes
- I. Micropipet and tips
- J. Centrifuge tubes or test tubes
- K. Purple Top Tubes

**22.5.2: Procedure** - The following standards are stored in the freezer and replaced as needed:Semen**A.** Preparation of semen standards

1. Aliquot 250µl volumes of neat semen into centrifuge tubes labeled with the sample-type, lot # (date of collection) and preparer's initials.
2. Store in a ziplock bag labeled with the sample-type, source (if available) and/or lot # (date of collection) and preparer's initials. (Note on the label if the sample is aspermic).

**B.** Christmas Tree and Sperm Hy-liter Stain Standards

1. Make a dilution (1:250 suggested) of neat spermic semen in dH<sub>2</sub>O and aliquot into centrifuge tubes labeled with the sample-type. Re-freeze remaining neat semen aliquot.
2. Store in a ziplock bag labeled with the sample-type, date/source of neat collection, lot # (date of preparation) and preparer's initials.

**C.** AP standards

1. Make a 1:10 dilution of neat semen in dH<sub>2</sub>O and aliquot 50µl volumes into centrifuge tubes labeled with the sample-type. Re-freeze remaining neat semen aliquot.
2. Store in a ziplock bag labeled with the sample-type, date/source of neat collection, lot # (date of preparation) and preparer's initials.

Urine (RSID™-Urine and p30 ABACard®)

- A. Saturate filter paper or swabs with neat female urine and dry overnight in the hood.
- B. Place sample made on filter paper into a coin envelope or repackage swabs into the paper sleeves and label with the sample-type, source (if available), lot # (date of preparation) and preparer's initials.

Fecal swabs (AP and Urobilinogen)

- A. Collect fecal material on swabs and dry overnight in hood.
- B. Re-package the swabs in the paper sleeves and label with the sample-type, source (if available), lot # (date of collection) and preparer's initials.

Breastmilk (p30 ABACard®)

- A. Store liquid breast milk in the freezer. Label with the sample-type, source (if available), lot # (date of collection) and preparer's initials.

- B. Thaw as needed and make a stain of the breast milk sample on filter paper. Re-freeze remaining sample.
- C. Dry overnight in hood and place into a coin envelope labeled with the sample-type, source (if available), date of collection, lot # (date of preparation) and preparer's initials.

Acid phosphatase standards

- A. Vaginal swabs
  - 1. Collect semen free vaginal samples (minimum of five days post coital) on swabs.
  - 2. Dry overnight in hood and re-package the swabs in the paper sleeves.
- B. Vaginal/semen mixed swabs
  - 1. Add 100µl of thawed 1:10 semen to each pre-made vaginal swab or collect post coital (~ 24 hours) vaginal/semen mixed swabs.
  - 2. Dry overnight in hood and re-package the swabs in the paper sleeves.
- C. Semen
  - 1. Make a stain on filter paper with 1:10 semen.
  - 2. Dry overnight in hood and place into a coin envelope.
- D. Oral swabs
  - 1. Collect oral sample on swabs.
  - 2. Dry overnight in hood and re-package the swabs in the paper sleeves.
- E. Urine Stain
  - 1. Saturate filter paper with urine.
  - 2. Dry overnight in hood and place into a coin envelope.
- F. Fecal swabs-described on the previous page
- G. Negative control-Use blank swabs, filter papers or cloth swatches as needed. Place into coin envelopes.
- H. Label each acid phosphatase standard with the sample-type, source (if available), lot # (date of collection/preparation) and preparer's initials. Store together in manila envelope labeled as 'AP' Standards.

Animal standards

- A. Collect blood samples from animal sources on white cloth swatches or filter paper and dry overnight in a designated area. Place into coin envelopes labeled with the sample-type, lot # (date of collection) and preparer's initials.

*Approved by Director: Dr. Guy Vallaro*

- B. Commercially available animal sera may be used as positive controls for the corresponding anti-sera. Aliquot 50µl volumes into centrifuge tubes labeled with the sample-type. Store in zip lock bags labeled with the sample-type, lot #, date received and preparer's initials.

**22.5.3: Procedure** - The following standards are maintained at room temperature and replaced annually (one (1) set of the expired standards are retained for research purposes and the remainder are discarded):

Blood (KM, o-Tolidine, Takayama, Ouchterlony, RSID™-Blood and HemaTrace®, crime scene kit)

- A. Collect blood in PTT's and make stains on filter papers. Refrigerate any remaining blood in the PTT labeled with the sample-type, source (if available), lot # (date of collection) and preparer's initials. Replace as needed.
- B. Dry overnight in a designated area.
- C. Cut and place each piece of stain into a coin envelope labeled with the sample-type, lot # (date of preparation), control date and preparer's initials.
- D. Replace the old standards with the new standards in the following Units: Forensic Biology, DNA, Trace (examiners and crime scene kit) and other examiners as necessary.

Semen (AP, p30 ABACard® and crime scene kit)

- A. Make a 1:10 dilution of thawed, neat semen in dH<sub>2</sub>O. Re-freeze remaining neat semen aliquot.
- B. Saturate each filter paper with approximately 1ml of the 1:10 dilution of semen.
- C. Dry overnight in hood.
- D. Cut and place each piece of stain into a coin envelope labeled with the sample-type, lot # (date of preparation), control date and preparer's initials.
- E. Replace the old standards with the new standards in the Forensic Biology Unit (examiners and crime scene kit).

Saliva (Phadebas®)

- A. Saturate filter papers with saliva.
- B. Dry overnight in hood.
- C. Cut and place each piece of stain into a coin envelope labeled with the sample-type, lot # (date of preparation), control date and preparer's initials.

- D. Replace the old standards with the new standards in the Forensic Biology Unit (examiners).

Blood, Semen, Saliva, Urine (alternate light sources)

- A. Make blood stains approximately 1" in diameter on black cloth swatches, ensuring that unstained substrate remains around the stain.
- B. Make separate saliva, urine and 1:10 semen stains approximately 1" in diameter on filter paper, ensuring that unstained substrate remains around each stain.
- C. Dry overnight in a designated area. Cut out the stains made on filter paper leaving unstained substrate around each.
- D. Check the new standard with the appropriate alternate light source(s) before use and record the results on the Body Fluid Standard Log Sheet.
- E. If appropriate for use, place each into a coin envelope labeled with the sample-type, lot # (date of preparation), control date and preparer's initials.
- F. If the appropriate results are not obtained, discard the standard, review the procedure and make a new standard.
- G. Replace the old standards with the new standards in the Forensic Biology Unit.

Negative controls

- A. Place blank filter paper into coin envelopes labeled with the sample-type, lot # (date of preparation), control date and preparer's initials.
- B. Replace the old standards with the new standards in the following Units: Forensic Biology, DNA and Trace (examiners).

**22.5.4: Procedure** - The following standards are stored at room temperature and replaced as needed:

Christmas Tree and Sperm Hy-liter control smears (if made in advance)

- A. Collect an epithelial cell (buccal) sample on a swab and form a smear onto a glass slide.
- B. With a micropipet, place approximately 3µl of thawed diluted spermic semen onto the smear. Re-freeze the remaining semen aliquot.
- C. Dry the positive control smear at room temperature or 37°C (do not apply open flame heat to the Sperm Hy-liter control smears).

- D. Label the smears with the sample-type, lot # (date of preparation) and preparer's initials and store in a slide box.

Sperm Hy-liter control swabs (made in advance as needed)

- A. Spermic semen/epithelial cell (buccal) swabs  
Collect epithelial cell (buccal) samples on swabs and add 25ul of thawed, neat semen (spermic) and dry overnight in hood. Re-freeze remaining semen aliquot.
- B. Separate spermic semen and epithelial cell (buccal) swabs  
1. Place 25ul of thawed, neat semen (spermic) onto swabs and dry overnight in hood. Re-freeze remaining semen aliquot.
2. Collect epithelial cell (buccal) samples on swabs and dry overnight in hood.
- C. Package the swabs in separate coin envelopes and label with the sample-type, lot # (date of preparation) and preparer's initials.

**22.5.5: Procedure** - The following blood standard is prepared for the LP Unit as needed:

- A. Label glass slides with the sample-type and lot # (date of preparation).
- B. Make fingerprint-type bloodstains from the PTT on the slides and dry overnight in a designated area.
- C. Forward to the LP Unit for use.

**22.6 REFERENCES**

GL-2 (Safety Manual)