

Document Title: Acid Phosphatase Reagent QC

Controlled: Yes, with red stamp present

Controlled By: Quality Manager

Prepared By: \_\_\_\_\_ Date: \_\_\_\_\_

Approved By: \_\_\_\_\_ Date: \_\_\_\_\_

**A. PURPOSE:**

To prepare reagents for semen screening test and to perform quality control on prepared reagents.

**B. RESPONSIBILITY:**

Forensic Science Examiners 1 and 2 in the Forensic Biology Section. Ordering information is maintained in a log book in the Forensic Biology Section.

**C. SAFETY:**

Use appropriate measures for the proper handling of acetic acid according to SOP-GL-2 (Safety Manual) and the Material Safety Data Sheets.

**D. PROCEDURE:**

1.  $\alpha$ -Naphthyl Phosphate Substrate Solution

a. Materials:

- aa. Acetate buffer (100ml):
  - Sodium acetate (crystal) 2g (\*1.23g if anhydrous)
  - Distilled water (dH<sub>2</sub>O) 100ml
- bb.  $\alpha$ -naphthyl phosphate (disodium salt) 0.187g
- cc. Glacial acetic acid
- dd. pH meter
- ee. Plastic tubes (12x75mm) and caps/parafilm
- ff. Test tube racks

b. Procedure:

- aa. Dissolve  $\alpha$ -naphthyl phosphate in 100ml of acetate buffer.
- bb. Uncover fill hole on electrode of pH meter.
- cc. Turn on pH meter.
- dd. Rinse electrode with dH<sub>2</sub>O and gently blot with kimwipe (don't rub or wipe).
- ee. Rinse electrode with sample prior to pH measurement.

- D. 1. b. ff. Place electrode into sample and wait for "pH" icon to stop flashing, read pH.
- ii. Bring to pH 5 by adding glacial acetic acid drop wise, allow to thoroughly mix and monitor the pH before adding another drop.
- jj. Aliquot 1.0 or 1.5ml volumes into plastic tubes in test tube racks and cover.

## 2. Fast Blue B Color Reagent

- a. Materials:
  - aa. Fast Blue B salt (o-dianisidine diazotate) 2.0g
  - bb. Distilled water (dH<sub>2</sub>O) 100ml
  - cc. Filter paper
  - dd. Plastic tubes (12x75mm) and caps
  - ee. Test tube racks
- b. Procedure:
  - aa. Dissolve Fast Blue B salt in dH<sub>2</sub>O.
  - bb. Filter if necessary.
  - cc. Aliquot 1.0 or 1.5 ml volumes into plastic tubes in test tube racks and cover.

## 3. Acid Phosphatase Reagent

- a. Test each new batch of reagent before use according to SOP-FB-11 (Screening Test for Semen). The results are recorded at the first indication of a pink or purple color change and observed for 15 seconds. Record the results according to the Acid Phosphatase Reagent Log Sheet.
- b. If the appropriate results are not obtained, discard the reagent, review the procedure and make new reagent.
- c. If the reagent is suitable for use, record the reagent, lot # (date of preparation) and examiner's initials on each rack and store in the freezer. AP reagent is suitable for use when non-semen containing body fluids that may contain acid phosphatase activity (i.e. vaginal secretions, fecal material and oral samples) yield weaker/slower results, or negative results, than semen containing body fluids (i.e. semen, semen/vaginal mixture). The PBS blank must yield a negative result. The use of AP reagent and the interpretation of AP results are addressed during training according to SOP-FB-11 (Screening Test for Semen).
- d. Discard any frozen aliquots after one (1) year.

**E. REFERENCES:**

1. Metropolitan Police Forensic Science Laboratory. Biology Methods Manual. 1978, pp.3-17 to 3-20.
2. SOP-GL-2 (Safety Manual).
3. Material Safety Data Sheets.