

Document Title: Acid Phosphatase Reagent QC

Controlled: Yes, with red stamp present

Controlled By: Quality Manager

Prepared By: Forensic Biology Section Date: _____

Approved By: _____ Date: _____

A. PURPOSE:

To prepare reagents for semen screening test and to perform quality control on prepared reagents.

B. RESPONSIBILITY:

Forensic Science Examiners 1 and 2 in the Forensic Biology Section. Ordering information is maintained in a log book in the Forensic Biology Section.

C. SAFETY:

Use appropriate measures for the proper handling of acetic acid according to SOP-GL-2 (Safety Manual) and the Material Safety Data Sheets.

D. PROCEDURE:

1. α -Naphthyl Phosphate Substrate Solution

a. Materials:

aa. Acetate buffer (100ml):

Sodium acetate (crystal)	2g (*1.23g if anhydrous)
Distilled water (dH ₂ O)	100ml

bb. α -naphthyl phosphate (disodium salt) 0.187g

cc. Glacial acetic acid

dd. pH meter

ee. Plastic tubes (12x75mm) and caps/parafilm

ff. Test tube racks

b. Procedure:

aa. Dissolve α -naphthyl phosphate in 100ml of acetate buffer.

bb. Uncover fill hole on electrode of pH meter.

cc. Turn on pH meter.

dd. Rinse electrode with dH₂O and gently blot with kimwipe (don't rub or wipe).

ee. Rinse electrode with sample prior to pH measurement.

- D. 1. b. ff. Place electrode into sample and wait for “pH” icon to stop flashing, read pH.
- ii. Bring to pH 5 by adding glacial acetic acid drop wise, allow to thoroughly mix and monitor the pH before adding another drop.
- jj. Aliquot 1.0 or 1.5ml volumes into plastic tubes in test tube racks and cover.

2. Fast Blue B Color Reagent

- a. Materials:
 - aa. Fast Blue B salt (o-dianisidine diazotate) 2.0g
 - bb. Distilled water (dH₂O) 100ml
 - cc. Filter paper
 - dd. Plastic tubes (12x75mm) and caps
 - ee. Test tube racks
- b. Procedure:
 - aa. Dissolve Fast Blue B salt in dH₂O.
 - bb. Filter if necessary.
 - cc. Aliquot 1.0 or 1.5 ml volumes into plastic tubes in test tube racks and cover.

3. Acid Phosphatase Reagent

- a. Test each new batch of reagent before use according to SOP-FB-11 (Screening Test for Semen) and the Acid Phosphatase Reagent Log Sheet. Record the required information.
- b. If the appropriate results are not obtained, discard the reagent, review the procedure and make new reagent.
- c. If the reagent is suitable for use, record the reagent, lot # (date of preparation) and examiner's initials on each rack and store in the freezer. **AP reagent is suitable for use when non-semen containing body fluids that contain acid phosphatase (i.e. vaginal secretions, fecal material and oral samples) yield weaker/slower results than semen containing body fluids (i.e. semen, semen/vaginal mixture). The PBS blank must yield a negative result. The use of AP reagent and the interpretation of AP results are addressed during training according to SOP-FB-11 (Screening Test for Semen).**
- d. Discard any frozen aliquots after one (1) year.

E. REFERENCES:

- 1. Metropolitan Police Forensic Science Laboratory. Biology Methods Manual. 1978, pp.3-17 to 3-20.

- E. 2. SOP-GL-2 (Safety Manual).
- 3. Material Safety Data Sheets.