

*Approved by Director: Dr. Guy Vallaro*

## **PART 1**

Turn on 56°C heat block

Make Tergazyme: 10mL dH<sub>2</sub>O + 0.5g Tergazyme (drawer) in 15mL Falcon tube

Clean grinder: Pour ~5mL Tergazyme into beaker, warm it on heat table (do not boil)

Clean grinder with swab/ warm Tergazyme, simulate grinding, rinse with dH<sub>2</sub>O

Repeat at least 2 times

Add 400μL 1 N H<sub>2</sub>SO<sub>4</sub> and let sit for 20 minutes (in acid hood)

Make SEB/DTT: 5mL SEB (pre-aliquoted) + 0.03g DTT (refrigerator) in 15mL Falcon tube

Stratalink: Falcon tube of- Tergazyme, Histoclear, 100% Ethanol, dH<sub>2</sub>O, SEB/DTT

10 screw-top tubes labeled- Tergazyme (2 tubes), Histoclear,

EtOH (2 tubes), dH<sub>2</sub>O (3 tubes), Q1, RBQ1

Stratalink for at least 15 minutes

Clean Dead Space Hood: Fresh 20% (in-house) bleach or 10% stabilized bleach, Kimwipe,

Falcon tube rack

Isopropyl wipe- tweezers, ruler, scissors/scalpel, 3 pipettors (p10, p200, p1000)

UV at least 15minutes

Dump out acid from grinder, rinse with dH<sub>2</sub>O

Pulse spin, pipette off excess liquid in hood

Stratalink grinder for at least 15 minutes

In dead space hood: Aliquot 1mL Terg., Histoclear, EtOH, and dH<sub>2</sub>O into labeled screw-top tubes

Add 200μL SEB/DTT to grinder, simulate grinding, transfer to RBQ1 tube

Add 1μL proK, place tube into 56°C heat block

Take out evidence and fill out Worksheets

*Approved by Director: Dr. Guy Vallaro*

In dead space hood: Measure hair in hood

Remove ~2cm from root end (if apparent) or more if needed (at the discretion of the analyst)

Transfer to 1mL Histoclear screw-top tube and sonicate for 20 minutes

Transfer hair to EtOH screw-top tube, invert 3-4 times

Transfer hair to dH<sub>2</sub>O screw-top tube, invert 3-4 times

Transfer hair to 1mL Tergazyme screw-top tube and sonicate for 20 minutes

Transfer hair to dH<sub>2</sub>O screw-top tube, invert 3-4 times

Transfer hair to new 1mL Tergazyme screw-top tube and sonicate for 20 minutes

Transfer hair to EtOH screw-top tube, invert 3-4 times

Transfer hair to dH<sub>2</sub>O screw-top tube, invert 3-4 times

Add 200µL SEB/DTT to grinder, transfer hair to grinder

Grind, pipette fluid to tube, add 1µL proK

Place into 56°C heat block for 2hrs up to Overnight

*Approved by Director: Dr. Guy Vallaro*

## **PART 2**

Stratalink: Rack with dH<sub>2</sub>O, 2 Microcon sets (filter and tube), 4 additional Microcon tubes

Clean Laminar Hood: Fresh 20% (in-house) bleach or 10% stabilized bleach, Kimwipe, Falcon tube rack

Isopropyl wipe- 2 pipettors (p200, p1000)

Autoclaved beaker to heat dH<sub>2</sub>O

UV at least 15 minutes

For the remainder of the procedure manipulate the RB tube and place into centrifuge before touching Q tube. Change gloves each time after handling Q tube. Apply UV to hood during spins.

Pulse spin tubes that have been incubating at 56°C to collect condensate

Add 200µL PCIA to each tube, vortex, and spin for 3 minutes at 10,000g (rcf)

Add 200µL of dH<sub>2</sub>O to Microcon set while waiting for tubes to spin

Pipette off supernatant and add to Microcon set, spin for 5 minutes at 3,000g (rcf)

Transfer filter to new Microcon tube, add 400µL of dH<sub>2</sub>O, and spin for 5 minutes at 3,000g(rcf)

Pour ~1mL dH<sub>2</sub>O into the beaker and heat until boiled (~3-5 minutes)

Add 60µL of the heated dH<sub>2</sub>O to the filter, invert into a new Microcon tube, vortex, and spin for 3 minutes at 10,000g (rcf)

Using a sterile pipette tip, determine the volumes of the RB and the sample extracts. The elution volumes shall be documented manually on QRM-4. The volume of the RB must not exceed the volume of the sample. If necessary, add dH<sub>2</sub>O to bring the sample up to the volume of the RB.

**mtDNA WI-03 Hair Extraction-Contaminated (blood) or  
mounted**

*Approved by Director: Dr. Guy Vallaro*

Document ID: 961

Revision: 1

Effective Date: 8/15/2014

Status: Retired

Page 4 of 4

RETIRED