

mtDNA WI-04 Hair Extraction-Contaminated (semen) or Consumption

Document ID: 962
Revision: 1
Effective Date: 8/15/2014
Status: Retired
Page 1 of 3

Approved by Director: Dr. Guy Vallaro

PART 1

Turn on 56°C heat block

Make Tergazyme: 10mL dH₂O + 0.5g Tergazyme (drawer) in 15mL Falcon tube

Clean grinder: Pour ~5mL Tergazyme into beaker, warm it on heat table (do not boil)
Clean grinder with swab/ warm Tergazyme, simulate grinding, rinse with dH₂O
Repeat at least 2 times
Add 400µL 1 N H₂SO₄ and let sit for 20 minutes (in acid hood)

Make SEB/DTT: 5mL SEB (pre-aliquoted) + 0.03g DTT (refrigerator) in 15mL Falcon tube

Stratalink: Falcon tube of- Tergazyme, Histoclear, 100% Ethanol, dH₂O, SEB/DTT
13 screw-top tubes labeled- Tergazyme (3 tubes), Histoclear, EtOH (3 tubes), dH₂O (4 tubes), Q1, RBQ1
Stratalink for at least 15 minutes

Clean Dead Space Hood: Fresh 20%(in-house) bleach or 10% stabilized bleach, Kimwipe,
Falcon tube rack
Isopropyl wipe- tweezers, ruler, scissors/scalpel, 3 pipettors (p10, p200, p1000)
UV at least 15 minutes

Dump out acid from grinder, rinse with dH₂O
Pulse spin, pipette off excess liquid in hood
Stratalink grinder for at least 15 minutes

In dead space hood: Aliquot 1mL Terg., Histoclear, EtOH, and dH₂O into labeled screw-top tubes
Add 200µL SEB/DTT to grinder, simulate grinding, transfer to RBQ1 tube
Add 1µL proK, place tube into 56°C heat block

Take out evidence and fill out Worksheets

In dead space hood: Place 1mL tube of Tergazyme into 56°C heat block for ~5-10 minutes

Measure hair in hood

mtDNA WI-04 Hair Extraction-Contaminated (semen) or Consumption

Document ID: 962

Revision: 1

Effective Date: 8/15/2014

Status: Retired

Page 2 of 3

Approved by Director: Dr. Guy Vallaro

Remove ~2cm from root end (if apparent) or more if needed (at the discretion of the analyst)

Transfer to 1mL of warmed Tergazyme tube, place tube back into 56°C heat block for 30 minutes

After the 30-minute incubation, transfer the hair to EtOH screw-top tube, invert 3-4 times

Transfer the hair to dH₂O screw-top tube, invert 3-4 times

Transfer hair to Histoclear screw-top tube and sonicate for 20 minutes

Place 1 of the additional tubes of Tergazyme into 56°C heat block to warm for ~5-10 minutes

Transfer the hair to EtOH screw-top tube, invert 3-4 times

Transfer the hair to dH₂O screw-top tube, invert 3-4 times

Transfer the hair to the new heated tube of Tergazyme and sonicate for 20 minutes

Place last additional tube of Tergazyme into 56°C heat block to warm for ~5-10 minutes

Transfer hair to dH₂O screw-top tube, invert 3-4 times

Transfer hair to remaining heated tube of Tergazyme and sonicate for 20 minutes

Transfer hair to EtOH screw-top tube, invert 3-4 times

Transfer hair to dH₂O screw-top tube, invert 3-4 times

Add 200µL SEB/DTT to grinder, transfer hair to grinder

Grind, pipette fluid to tube, add 1µL proK

Place into 56°C heat block for 2hrs up to Overnight

mtDNA WI-04 Hair Extraction-Contaminated (semen) or Consumption

Document ID: 962
Revision: 1
Effective Date: 8/15/2014
Status: Retired
Page 3 of 3

Approved by Director: Dr. Guy Vallaro

PART 2

Stratalink: Rack with dH₂O, 2 Microcon sets (filter and tube), 4 additional Microcon tubes

Clean Laminar Hood: Fresh 20%(in-house) bleach or 10% stabilized bleach, Kimwipe, Falcon tube rack

Isopropyl wipe- 2 pipettors (p200, p1000)

Autoclaved beaker to heat ddH₂O

UV at least 15 minutes

For the remainder of the procedure manipulate the RB tube and place into centrifuge before touching Q tube. Change gloves each time after handling Q tube. Apply UV to hood during spins.

Pulse spin tubes that have been incubating at 56°C to collect condensate

Add 200µL PCIA to each tube, vortex, and spin for 3 minutes at 10,000g (rcf)

Add 200µL of dH₂O to Microcon set while waiting for tubes to spin

Pipette off supernatant and add to Microcon set, spin for 5 minutes at 3,000g (rcf)

Transfer filter to new Microcon tube, add 400µL of dH₂O, and spin for 5 minutes at 3,000g(rcf)

Pour ~1mL dH₂O into the beaker and heat until boiled (~3-5 minutes)

Add 60µL of the heated dH₂O to the filter, invert into a new Microcon tube, vortex, and spin for 3 minutes at 10,000g (rcf)

Using a sterile pipette tip, determine the volumes of the RB and the sample extracts. The elution volumes shall be documented manually on QRM-4. The volume of the RB must not exceed the volume of the sample. If necessary, add dH₂O to bring the sample up to the volume of the RB.