

*Approved by Director: Dr. Guy Vallaro***Amplification – Examiner Setup**

The ‘**Examiner Setup**’ module is used to make amplification decisions and will be completed by an authorized DNA analyst.

1. Open **Processing → DNA Processing → Examiner Setup**
2. The **Examiner Setup** module can be filtered for ease of view. Suggested viewing is as follows:
  - a. Right click → Select Show Group By. Drag the Batch Barcode column to the top. This will separate samples by Extraction batch ID.
  - b. Right click → Select Filter Type header. Use this feature to select one detector to view at a time (particularly step 4).
  - c. Right click → Select Save Layout to maintain this view for all future viewings.
3. Isolation batches available for amplification will be displayed in the **Quantitation Results** worklist.
4. Each sample will be displayed 3 times – once for each detector value. However, only one detector value needs to be selected for subsequent testing. The detector selected will be the quantification value used for the amplification calculations. RBs that were not quanted will be listed only once with no quant value.

The Action column recommends “use as is” vs. “dilution” based on target amplification values


Bar Code ▾	Lab Case Num ▾	Exhibit Number ▾	Lot Number ▾	Processing Analyst ▾	Report Analyst ▾	Detector ▾	Avg Auto ▾	Quantity (ng/ul) ▾	Action ▾	Extract Remaining (ul) ▾	Exhibit Description
SP016.A	DSS-25-000010	001-003-02		Julia Gordon		T.Large Autosomal	1.78415	2.8493	Dilute	27.0000	Vaginal swabs
SP016.A	DSS-25-000010	001-003-02		Julia Gordon		T.Small Autosomal	1.78415	0.719	Dilute	27.0000	Vaginal swabs
SP016.A	DSS-25-000010	001-003-02		Julia Gordon		T.Y		1.2532	Dilute	27.0000	Vaginal swabs

5. Samples highlighted in crimson (IPC=Yes) are an indication of a possible inhibited sample.
6. Samples highlighted brown are samples with a Degradation Index  $\geq 3$ . Select the large autosomal quant values for these samples.
7. If necessary, any samples may be reworked by selecting the sample and clicking **Rework**. Select the appropriate entry point and add comments.

*Approved by Director: Dr. Guy Vallaro*

Well	Lab Number	Exhibit Number	Bar Code	Lot Number	Processing Analyst	Report Analyst
A02	DSS-25-000010	001-003-02	SP016.A		Julia Gordon	

8. To halt a sample, select the sample (any detector line) and click the **Discontinue** button. A discontinue reason will be added using the autofill option (Right click → Autofill → Examiner Setup comments).



9. Adding a sample to an amplification kit: While selecting the appropriate detector for calculations, highlight the samples for amplification. (ex: highlight the large autosomal value for degraded samples, Y value for Y-STR testing, etc. )
10. Select the amplification kit from the kit dropdown and click **Add**. Do not forget to add the appropriate blanks and EP1 as well.

*Approved by Director: Dr. Guy Vallaro*

Kit: Globalfiler  
Globalfiler Express  
Yfiler Plus

Add Quantity: 0.0000

Number	Bar Code	Lot Number	Extract Remaining
--------	----------	------------	-------------------

11. Adjust the **Action** column, **Dilution** and/or **Reaction DNA** columns as needed. Calculations will be set to target 15uL of input DNA volume. Leaving the target dilution concentration at 0.02 ng/uL. This can be changed to 0.3ng/ul to decrease the volume of total dilution if necessary.

Quantity: 0.0200

Amplification Setup Sheet : 1

Lab Case Number	Exhibit Number	Bar Code	Action	Quantity (ng/uL)	Kit	Dil. DNA (uL)	Dil. Buffer (uL)	Dil. Conc (ng/uL)	Reaction DNA (uL)	Reaction Diluent (uL)	DNA Amplified (ng)	Ex
DSS-25-000010	001-003-02	SP016.A	Dilute	0.719	Globalfiler	1	35	0.02	15	0	0.30	27

12. Samples may be added to the Amplification Setup Sheet window multiple times (ex: amplifying a sample in both GF and YFP) by selecting different detector values.

**Note: If a sample is to undergo multiple amplifications, this must be done prior to marking the sheet complete.**

Select a sample choose the other kit and click add.

13. Click **Mark Complete**.

### Amplification – Setup & Processing

1. Open **Processing** → **DNA Processing** → **Amplification Setup**.
2. Click **Create Batch**. The **Amplification Batch Create** screen opens (another window).
3. Click **Create**. Select a **Batch Type**. The Fill order is not relevant to manual processing.
4. Click **Create** and scan the generated amplification batch bar code.
5. Select the **Source Batch** and **Allocate** to the **Destination Amplification Batch**. Samples can be moved around after hitting save. Be sure if moving the negative control, you choose “consumable negative control” and not “amplification blank”.
6. Click **Complete** when all samples are added.

*Approved by Director: Dr. Guy Vallaro*

7. In the **Amplification Setup** module which is already open, select/scan the batch to be processed and click **Select Scenario**.
8. Scan in all consumables, batch barcodes and hood/workstation. The consumables window will automatically calculate the reagent volumes needed to make the master mix.

Consumables : 1						
Group Type	Group Name	Extra	Extra Number			
All		# Extra Reactions	2			
Optional	Fixed	Consumable Name 1	Format 1	Bar Code 1	Reagent Amount	Measurement
<input type="checkbox"/>	<input type="checkbox"/>	GlobalFiler MasterMix	GM#####	GM000005	52.50	ul
<input type="checkbox"/>	<input type="checkbox"/>	GlobalFiler Primer Mix	GP#####	GP000005	17.50	ul
<input type="checkbox"/>	<input checked="" type="checkbox"/>	dH2O	WA#####	WA000001		
<input type="checkbox"/>	<input checked="" type="checkbox"/>	DNA Control 007 (0.1ng/μL)	GC#####	GC000003	3.00	ul

9. For ease of aliquoting, volumes and dilutions can be viewed by clicking on the Content button on the right-hand column. Select the **Amp Setup Sheet Results** tab and organize columns as needed. This window can remain open throughout the process, separate from the Amplification Setup window.

Batch Content (BatchDetailForm)

Batch ID: 5557    Bar Code: GFA-250312-59    Status: Amplification Batch Create

Batch Content    Quant Results    Amp Setup Sheet Results

Amplification Results : 8

Bar Code	Lab Case N	Exhibit Num	Dil. DNA (uL)	Dil. Buffer ( ▾ )	Quantity (ng/uL)	Dil. Quantity (ng/uL)	Reaction DNA (uL)	Reaction Buffer (uL)	DNA Amplified (ng)
CC010.B	CLC-25031...	1-2	1.0	370.0	7.4028	0.02	15.0	0.0	0.30
CC009.B	CLC-25031...	1-1	1.0	254.0	5.1051	0.02	15.0	0.0	0.30
CC012.A	CLC-25031...	1-2	1.0	235.0	4.7305	0.02	15.0	0.0	0.30
CC011.A	CLC-25031...	1-1	1.0	124.0	2.5007	0.02	15.0	0.0	0.30
250033.A					0.0		15.0	0.0	0.00
250034.A					0.0		15.0	0.0	0.00
250033.B					0.0		15.0	0.0	0.00
250034.B					0.0		15.0	0.0	0.00

10. When ready to aliquot, click **Start Process**.
11. Complete the tube check by scanning the tube labels.

Note: Scan the DNA Control 007 and dH2O lot numbers as the positive and negative, respectively.

*Approved by Director: Dr. Guy Vallaro*

12. When the amplification setup has been completed, click **Complete Process** and record the results using the Complete Batch Activity screen. Click **Save**.
  - **Process Successful:** the batch advances to the next processing step.
  - **Process Aborted:** the batch remains on the Extraction Batches worklist. A **Batch Comment** is required with this option.
  - **Process Failed:** the batch is abandoned and all samples return to Quant Setup. A batch comment is required with this option.
13. Document the thermocycler used. Open **Processing → DNA Processing → Amplification**. This can be done in the post-amp lab.
14. Scan the batch from the Batches Worklist.
15. Scan the thermocycler bar code.

Note: if a C-1000 is used, indicate in the batch comments which block is used.
16. Click **Start Process**.
17. Fill out the **Complete Batch Activity** window by completing the process when done.
18. Next step is STACS SOP-7 CE.