

BATCH PAPERWORK & ANALYSIS**Analyst**

- 1) Bring all data files into GeneMarker and analyze using GF_CTDSS.
- 2) Save GeneMarker project to file folder that HID files are in (U:\3130 NUCLEAR DNA CASEWORK & QC\To Be Analyzed).
- 3) Analyze samples, adding comments for all click-offs. Click-offs should mostly be dye-artifacts, raised baseline, and pullup peaks.
- 4) If a sample needs to be re-injected, or if for another reason you are not reporting out a sample, disable the sample by right-clicking on the file name in the left panel and note the reason on QR-325
- 5) Re-save GeneMarker file when all click-offs have been made and all samples have been disabled.
- 6) Print out e-grams from samples and controls.
- 7) Export report as a text file to your "Analysis Export File" in your "To Be Analyzed" folder (this text file is used by the Concordance Checker, Cross-Comparison/Staff Search, Contributor Estimator, Deconvolution workbook, and Allele Table macros).
- 8) Run concordance checker (DNA QR-37) if positive controls are present, using text file created.
- 9) Repeat steps 1 thru 8 if there are multiple injection sheets.
- 10) Complete batch cross-comparison.
- 11) Initial and date QR-4A.
- 12) Give batch to Technical Reviewer.

Technical Reviewer

- 1) Review all batch paperwork.
- 2) Open the analyst's fully analyzed project.
- 3) Review samples, note any differences to the 1st Analyst. Check that:
 - a. Correct analysis template is used.
 - b. Analyst's assessment of size standards, ladders, controls, and potential contamination are appropriate.
 - c. All called peaks appear to be true alleles.
 - d. All deletions/edits made are appropriate.
- 4) Initial and date QR-4A.
- 5) Give batch back to analyst.

Analyst

- 1) Make changes to project(s), if necessary, according to the technical review.
- 2) For each sample project, start a new project by bringing in the same data files into GeneMarker, and analyze them using GF_CTDSS_STRmix template. Save GeneMarker project with identical name as before, but with "STRmix" added to the end.

Approved by Director: Dr. Guy Vallaro

- **OR** you can re-analyze your data from the current project by pressing the green play button again. A message will come up asking if you want to re-call the size standard. You want to check the box for “Call Size Again” and click on “Apply to all”. Don’t forget when you save this file to save it with STRmix added to the end.
- 3) Using your analysis on the e-gram peak tables as a guide, click off all the same peaks you just clicked off with the following exceptions (that will all come up as non-concordances when you run your project comparison, to be noted in the comments section):
 - a. If you clicked off a stutter peak because pullup contributed to it, leave this peak in your STRmix project, provided that the peak is still greater than AT, with pullup subtracted, and the pullup is not responsible for the majority of the peak height.
 - b. If you clicked off a peak as high stutter, leave this peak in your STRmix project.
 - c. Same goes if any of these peaks are a combo of different stutters
 - i. There might be exceptions to this, for example, if the subtraction of pull-up brings the stutter peak below AT, STRmix should not see this peak, and if the pull-up contribution adds an unusually large amount of height to the stutter peak, it might be better for STRmix not to see such a high peak. Ask for assistance if you are unsure.
- 4) Disable all samples that you disabled in the previous project.
- 5) Re-save GeneMarker file when all edits have been made and all samples have been disabled.
- 6) Export report as a text file, same location as file created in step 7 on page 1.
- 7) Run DNA-QR-301 Project Comparison Tool using the exported text files for your fully analyzed project and your STRmix project.
 - a. If there are any discrepancies not accounted for in your analysis, you must correct them in the appropriate GeneMarker project. Then re-save the GeneMarker projects AND re-export the text files AND re-run project comparison macro.
 - b. Keep repeating this until the number of discrepancies can be explained by your analysis.
 - c. Print out this sheet for the batch paperwork; you only need the final copy in your batch.
- 8) Create DNA-QR-302 Contributor Estimation Worksheet for each sample with results. If it’s a single source sample, or you feel you don’t need the macro to assist you in determining number of contributors, you still need this sheet, but add notes to the comments area as to your determination of number of contributors.
 - a. Note minimum number of male contributors on designated line.
- 9) Initial and date QR-4A (Batch Paperwork Review Worksheet).
- 10) Give batch to Technical Reviewer to review additional worksheets.

Technical Reviewer

- 1) Review all additional paperwork (Project Comparison Worksheet, Contributor Estimation Worksheets).
- 2) Initial and date QR-4A.
- 3) Give to Administrative Reviewer

*Approved by Director: Dr. Guy Vallaro***Administrative Reviewer**

- 1) Review all batch paperwork.
- 2) Initial and date QR-4A.
- 3) Give batch back to Analyst.

Analyst

- 1) Review & make any corrections.
- 2) Initial and date QR-4A. Return batch to Technical Reviewer.

Technical Reviewer

- 1) Review any changes made by 1st analyst.
- 2) Add batch name to spreadsheet on U: Drive.
- 3) Initial and date QR-4A.
- 4) Return batch to Analyst.

Analyst

- 1) Make photocopies of injection sheets, Quant Trio Report, and halt at quant sheets for appropriate case jackets in batch. Original sample e-grams & original contributor estimator worksheets are just for case jackets, not to be in batch paperwork.
- 2) Scan batch paperwork (not containing sample e-grams & contributor estimator sheets) & save PDF file to U: Drive.
- 3) Move HID and project files from "To Be Analyzed" to "Completed" folder on U: Drive.
- 4) Give batch paperwork to Technical Reviewer to review.

Technical Reviewer

- 1) Ensure all pages scanned and PDF file on U: Drive is complete and accurate.
- 2) Ensure HID folders have been moved to U:\3130 NUCLEAR CASEWORK & QC\Completed.
- 3) Add initials to batch paperwork spreadsheet.

CASEWORK**Analyst**

- 1) Enter all possible profiles into CODIS.
- 2) Make comparisons of all knowns to all Q's in case.
 - a. Does a sample have 5+ contributors?
 - i. If so, do any of the knowns CBE to the sample warrant being conditioned? If yes, does that bring it down to 4 unconditioned contributors?
 1. Yes – can interpret.

Approved by Director: Dr. Guy Vallaro

2. No – Clear major present – can make comparison to deduced major profile. No comparisons will be made to minor components of profile.
3. No – No clear major present – No comparisons will be made to profile.
- b. Samples that are clear eliminations are not run through STRmix, “Eliminate by Analyst” in report macro.
- c. Comparisons that are not clear eliminations, are run through STRmix.
 - i. Reported out comparison based on LR calculated.
 - ii. The “Sample ID” field shall contain both item # and known # (if applicable), so that there is a unique identifier in the footer of all pages when printed.
- d. If reporting out a statistic for a comparison, deconvolute in STRmix with appropriate knowns conditioned to profile.
 - i. Report out comparison based on LR calculated.
 - ii. If not reporting out a statistic (i.e. intimate sample with no suspect known submitted, single source sample when known expected on sample), no need to run STRmix, enter “assumed” in report workbook. STRmix may be run if necessary to determine if assuming is appropriate.
- e. Does a sample only have one allele at one locus? → inconclusive or elimination by analyst
- f. Are you unable to determine the number of contributors with reasonable scientific certainty?
 - i. Clear major present – Can make comparison to deduced major profile. No comparisons will be made to minor components of profile.
 - ii. No clear major present – No comparisons will be made to profile.
- g. Is there a problem with the deconvolution discovered Re: 2nd diagnostics that can’t be resolved? → inconclusive or elimination by analyst
- 3) If a sample is deconvoluted through STRmix:
 - a. The file will automatically save onto the server, under F:\results.
 - b. Review the summary report. Ensure that the input file is correct. The letter “S” has previously been added to the beginning of the sample names of the STRmix input file when the project comparison macro was run.
 - c. Print out pages 1 thru 3 of the STRmix report.
 - d. Scrutinize summary reports for any problems with the deconvolution.
 - i. Fill out DNA-QR-303 for the 2nd diagnostics for each deconvolution.
 - e. Move the file in the F:\results folder to your folder in F:\results.
 - f. Create a folder in your folder for this case, named with the case number in format DSS-XX-XXXXXX, and move that folder into this newly created case specific folder.
- 4) Review case jacket, write report.
- 5) Give case jacket to TR.

Technical Reviewer

Approved by Director: Dr. Guy Vallaro

- 1) Review case jacket, comparisons and report as you normally would.
 - a. If you disagree with an eliminated comparison, the 1st Analyst will have to run this sample through STRmix.
 - b. Review the summary report. Ensure that the input file is correct. The letter “S” has previously been added to the beginning of the sample names of the STRmix input file when the project comparison macro was run.
- 2) If you need to review more of the summary report than is in the case jacket, you can find it on the STRmix server, F:\Results in the Analyst’s folder, in the folder specific to this case.
- 3) Give case jacket to AR.

Administrative Reviewer

- 1) Review case jacket and report as you normally would.
- 2) Give case jacket back to Analyst.

Analyst

- 1) Make all appropriate corrections to report and case jacket.
- 2) If the TR disagreed with one of your elimination calls, you must run that sample through STRmix, print out appropriate pages and fill out DNA-QR-303.
 - a. The report will have to be re-worked:
 - i. If it is still an elimination, it is now “elimination by STRmix”, and not “elimination by analyst”, report wording will be different.
 1. If it is no longer an elimination, report out result based on LR.
- 3) Print out final copy of report.
- 4) If there are STRmix results, move the folder for this case from your folder on the server’s F: Drive to F:\results\Completed.
- 5) Give case jacket to TR for correction check.

Technical Reviewer

- 1) Ensure all corrections are accurate.
- 2) If additional deconvolution work needed to be completed, make sure new wording in report is accurate and review all additional paperwork.
- 3) If there’s STRmix data, check to make sure the folder has been moved to F:\results\Completed.
- 4) Give case jacket to AR for correction check.

Administrative Reviewer

- 1) Ensure all corrections are accurate.
- 2) Give case jacket back to analyst.

CODIS

*Approved by Director: Dr. Guy Vallaro***Analyst**

- 1) If there is a CODIS hit for a sample of yours, you will be notified by the CODIS administrator/alternate or their designee.
- 2) The CODIS hit might warrant the sample to be deconvoluted in STRmix. If a deconvolution of this sample has already been completed, the case jacket can be brought to the CODIS administrator/alternate or their designee.
- 3) If a deconvolution of this sample has not yet been completed, deconvolute the sample, making sure to condition any knowns that you assumed when writing the report.
 - a. Print out page 1 of the STRmix report, fill out DNA-QR-303, and give the case jacket to your TR to review the additional paperwork. You can check, initial, and date next to 3.D. on DNA-QR-4.

Technical Reviewer

- 1) Review additional paperwork in case jacket.
- 2) Check, initial and date next to 3.D. on DNA-QR-4.
- 3) Give the case jacket to the AR

Administrative Reviewer

- 1) Review additional paperwork in case jacket.
- 2) Document review on DNA-QR-4.
- 3) Give case jacket to the CODIS administrator/alternate or their designee.

Analyst

- 1) When you receive your case jacket back for initials, move your STRmix folder to F:\results\completed. Check, initial, and date that you did so on DNA-QR-4.

OTHER

- 1) If a known comes in later on that warrants being conditioned, all deconvolutions need to be repeated with a new deconvolution, and a new statistic is calculated.
- 2) If a known comes in later that's not clearly eliminated regarding a sample that has already been deconvoluted, use the "LR to previously analysis" function of STRmix, and do not repeat the deconvolution.
- 3) Similarly, if you have a Q in the case, and two knowns are CBE to that sample, only deconvolute the sample once. For the 2nd known, use the "LR to previous analysis" function of STRmix.

Approved by Director: Dr. Guy Vallaro