GL 21.1 Balance Performance Check

CT DESPP Division of Scientific Services

Revision: 3

Page 1 of 1

Dolongo	Performance	Chaol	1-
Dalance	Performance	Chec	Κ

Balance	ID·	
Darance	112.	

Effective: 1/22/2025

	Low:	Mid:	High:	Other if Required:		
Reference Weight ID						
Replicate						
1						
2						
3						
4						
5						
6						
7						
8						
9						
10						
Average						
%Accuracy						
Standard						
Deviation Coefficient of						
Variation						
Analyst:		Date:	Unit:			
☐ Based on the performance of the acceptability criteria above the balance meets the criteria						
for use.						
☐ Based on the performance of the acceptability criteria above the balance does not meet						
the criteria for use. The balance has been marked 'do not use'.						
,						

Acceptability criteria (for balances with a max load of \leq 700g): %Accuracy +/-2%, SD <1, Coefficient of Variation \leq 2% * Acceptability criteria (for balances with a max load of \geq 701g): %Accuracy +/-5%, SD <1, Coefficient of Variation \leq 3% * *Units may require tighter tolerances based on the typical use of the balance.

- 1. Percent Accuracy: (expected value average value) (expected value)*100
- 2. Standard Deviation:

$$\sigma = \sqrt{\mathrm{E}((X - \mathrm{E}(X))^2)} = \sqrt{\mathrm{E}(X^2) - (\mathrm{E}(X))^2}$$

3. Coefficient of Variation: $V = (\sigma / \text{average value})*100$

If an excel spreadsheet is used to calculate the %Accuracy, SD or Coefficient of Variation attach a copy of the spreadsheet for the reviewer.