## Weight Performance Check

vi eight i ci ioi munce en	CCIX				
Expected Weight					
Weight Serial Number or					
Identifier					
"As left" value from last external					
calibration certificate					
Traceable Report Number					
Balance used					
Replicate 1					
2					
3					
4					
5					
6					
7					
8					
9					
10					
Replicate Average					
% Accuracy <sup>1</sup>					
Standard Deviation					
Coefficient of Variation					
*Acceptable for Use?					
Analyst:	Date:				
Reviewed by:	Date:_				
-					
NT		, 11 0	111.1	. 1 1	. 1
Note: Any weight which is determined to be unacceptable for use will be immediately taken out of					
service and reported to the Unit Supervisor and or Quality Section.					

Revision: 4

\*Acceptability criteria % Accuracy (+/-2%) SD, <1, Coefficient of Variation <3%,

Class 1 or Ultraclass weights  $\leq$ 2% of expected weight as per the last calibration certificate other classes  $\leq$  5%

1. Percent Accuracy: (expected value – average value)
(expected value)\*100

<sup>1</sup> For expected value use "as left" value from certificate

2. Standard Deviation:

$$\sigma = \sqrt{E((X - E(X))^2)} = \sqrt{E(X^2) - (E(X))^2}$$

3. Coefficient of Variation:  $V = (\sigma / \text{ average value})*100$ 

If an excel spreadsheet is used to calculate the %Accuracy, SD or Coefficient of Variation attach a copy of the spreadsheet for the reviewer.

Maintain this record in the Unit Logbook (paper or electronic) for this device.