

## Vibrio Isolate or Culture

Revised: 1-31-2025

<b>Test Description</b>	Confirmatory identification of enteric bacteria of public health significance, <i>Virbio species</i>
<b>Test Use</b>	Identification and serotyping of <i>Virbio species</i>
<b>Test Department</b>	Bacteriology Phone: (860) 920-6596 FAX: (860) 920-6721
<b>Methodology</b>	Bacterial culture on selective media, biochemical identification, API 20E, serotyping
<b>Availability</b>	Daily, Monday - Friday
<b>Specimen Requirements</b>	<ul style="list-style-type: none"> <li>• Pure culture of enteric pathogen on agar slant (preferred for shipping) or plate.</li> <li>• Stools tested positive for <i>Vibrio species</i></li> <li>• Freshly passed stool specimen collected early in the course of disease</li> </ul>
<b>Collection Kit/Container</b>	To obtain collection kit, refer to Collection Kit ordering Information.
<b>Collection Instructions</b>	<p>Stool specimens from patients with diarrheal illness must be collected as early in the course of the disease as possible. Collect into a clean, dry container.</p> <p>Transfer specimen to Cary Blair transport media. Fill to indicated line (15 mL).</p>
<b>Specimen Handling &amp; Transport</b>	Transport isolates to the laboratory at room temperature (15-25°C). Avoid temperature extremes. Follow all applicable federal packaging and shipping regulations. Stool specimen must arrive within 7 days of collection
<b>Requisition Form</b>	Clinical test requisition OL-9B
<b>Required Information</b>	<ul style="list-style-type: none"> <li>• Name and address of submitter. Two patient identifiers (ie. name, DOB, Acc.#, MRN), town of residence (city, state, zip). Specimen type or source, date collected, test requested, identity of the enteric pathogen to be confirmed</li> <li>• Ensure information on the requisition matches the specimen.</li> </ul>
<b>Unacceptable Conditions</b>	<ul style="list-style-type: none"> <li>• Unlabeled specimens</li> <li>• Specimen with less than 2 unique identifiers (ie. name, DOB, Acc.#, MRN)</li> <li>• Specimens that leak or containers that have broken in transit</li> <li>• Stool specimen received &gt;7 days of collection</li> <li>• Non-viable isolates</li> </ul>