

MYCOBACTERIOLOGY TESTING SERVICES GUIDE

**DR. KATHERINE A. KELLEY
CONNECTICUT DEPARTMENT OF
PUBLIC HEALTH LABORATORY**



Connecticut Department of Public Health
Mycobacteriology Laboratory

395 West Street, Rocky Hill, CT 06067
(860) 920-6649

Introduction

The Dr. Katherine A. Kelley Connecticut Department of Public Health Mycobacteriology Laboratory serves as a resource for hospitals and health clinics in confirming and identifying *Mycobacterium tuberculosis* complex (MTBc) and other mycobacteria of clinical interest. The Mycobacteriology Laboratory receives and processes several hundred specimens a year while providing accurate results in a short time. The quick turnaround of results ensures that submitters and clinicians are prepared to provide the most appropriate treatment for their patients. This guide will outline the process of how to submit a specimen for testing, as well as serve as an educational tool in the procedures conducted to confirm that a specimen harbors MTBc.

Who We Are

Bacteriology/Bioterrorism Response/Mycobacteriology Department



Rear, Left to Right : Dave Santoro, Hongli Dong, Mary Anne Banevicius, Beata Harvey, Doneisha Rumble, Rik Martinez
Front, Left to Right : Mark Harkins, Jocelyn Florano Boguszewski, Christina Nishimura, Barbara Wang

Not Pictured : Nick Brunetti

Connecticut Department of Public Health Laboratory Director

Jafar Razeq, PhD, HCLD (ABB)

Infectious Diseases Division Director

Anthony Muyombwe, PhD, HCLD (ABB)

Bacteriology/Bioterrorism Response/Mycobacteriology Supervisor

Christina Nishimura

Specimen Submission Procedures

Timely and accurate testing is contingent upon proper specimen submission procedures being exercised by the submitter or clinician. This includes appropriate **specimen collection** and **shipping**, in addition to a properly filled **clinical test requisition form**.

Specimen Collection

Correct specimen collection is a critical first step in the analytic process. Appropriate shipping conditions also ensure specimen quality is maintained throughout transport to the Laboratory.


Acceptable Specimens and Shipping Conditions

Specimen Type	Shipping Container	Acceptable Testing Volume	Shipping Conditions
Pulmonary Specimens (sputum, bronchial wash, bronchial lavage)	TC OL40 or similar container	2-5 mL	Ice Pack
Extra-Pulmonary Body Fluids	TC OL40 or similar container	2-5 mL	Ice Pack
Tissues and CSF	TC OL40 or similar container	Tissue samples, biopsies; CSF recommended minimum volume 0.25 mL	Ambient temperature
Blood and Bone Marrow	Blood collection tubes containing heparin (green top) or sodium polyanethol sulfonate SPS (yellow top)	2-5 mL	Ambient temperature
Culture Isolates (Referred Cultures)	Lowenstein-Jensen agar or Middlebrook 7H10 and 7H11 agar; liquid media from automated test systems (i.e., BACTEC MGIT broth)	2-5 mL or solid media	Ambient temperature

PLEASE NOTE:

- Swabs are strongly discouraged as a clinical specimen source.
- The shipment of plate cultures is discouraged.
- Insufficient volume of specimen may not allow for complete testing.
- Extreme overgrowth of non-acid fast organisms may render a specimen unusable.

Unacceptable Specimens:



- Unlabeled or illegibly labeled specimens
- Broken or leaking specimens
- Blood and bone marrow collected in tubes other than heparin or SPS
- Specimens without two distinct patient identifiers

Clinical Test Requisition

A **clinical test requisition** (OL-9B) must accompany every specimen (Appendix A) and is also available electronically – copy and paste this link into your browser:

https://portal.ct.gov/-/media/Departments-and-Agencies/DPH/laboratory/labhome/lab-forms/ClinTestReq_OL9B_FILL-Up1-10-24.pdf

The following information is required on the clinical test requisition (as indicated in Appendix A):

- Name and address of submitter (and/or profile number)
- Patient name or unique identifier, date of birth and address
- Specimen type or source of collection
- Date of specimen collection
- Testing requested
- This information **must** match the labeling on the specimen itself. **Two distinct patient identifiers** that match between the specimen labeling and the requisition information are required.
- Any missing or incorrect information may delay the release of results until the issue is resolved.
- Any specimens that are unlabeled or illegibly labeled will be rejected for testing.

Cepheid GeneXpert Nucleic Acid Amplification Testing (NAAT)

If **nucleic acid amplification (NAA) testing** is requested (please see special specimen requirements regarding this request below), **an additional form** must be filled out by the submitter or clinician (Appendix B). This form is available electronically – copy and paste this link into your browser:

<https://portal.ct.gov/-/media/Departments-and-Agencies/DPH/laboratory/labhome/lab-forms/MTB-NAAT-REQUEST-FORM-CTDPH-LAB-0418.pdf>



WHAT IS NAA TESTING AND ARE THERE SPECIAL SPECIMEN REQUIREMENTS?

Nucleic acid amplification (NAA) testing (using the Cepheid GeneXpert) is used to identify the presence of MTBc and potential resistance to rifampin, an important first-line drug. This test is offered for **unprocessed** sputum, bronchial wash and bronchial lavage **clinical specimens only**. This testing requires that the **patient has received less than three days of antituberculous therapy**, and that the **specimen is collected less than ten days before receipt in the laboratory**.

Collection Kits and Shipping Materials

Specific **specimen collection kits** for mycobacterial clinical specimens (TC OL40) and shipping materials are available free-of-charge by the Department of Public Health Laboratory Outfit Room.

For collection kits and shipping supplies, contact the DPH Laboratory Outfit Room at: (860) 920-6674 or (860) 920-6675, or by email at dph.outfitroom@ct.gov

ADDITIONAL SHIPPING INFORMATION

Shipping of cultures known or suspected of containing *Mycobacterium tuberculosis* complex must be packaged and shipped in compliance with “Category A Infectious Substances” guidelines by appropriately trained personnel. Please refer to the following for additional guidance:

<https://www.phmsa.dot.gov/sites/phmsa.dot.gov/files/2020-04/Transporting-Infectious-Substances-Safely.pdf>

<https://shop.saftpak.inmarkinc.com/products/shipping-category-a-biological-substance-and-related-materials-training-course>

Mycobacteriology Lab Services Available

METHOD	PURPOSE	TIMEFRAME FOR RESULTS (From Date of Receipt)
Acid-fast bacilli (AFB) Smear	<p>Determine the presence or absence of mycobacteria in a clinical specimen.</p> <p>Monitor response to drug treatment.</p> <p>Characterize acid fast organisms in a culture isolate.</p>	<p>Within 24 hours (except for specimens received Friday after noon)</p> <p>Positive clinical smears from a new patient are communicated to the submitter as <i>critical values*</i>.</p>
Nucleic Acid Amplification Testing (NAAT)	<p>Determine the presence of MTBc DNA in a clinical pulmonary specimen, and the presence or absence of genetic markers for rifampin resistance, a potential indicator of multi-drug resistant (MDR) MTBc.</p>	<p>Within 24-48 hours (except for specimens received Friday after noon)</p> <p>All results, positive or negative, are communicated to submitters as <i>critical values*</i>.</p>
Culture and Identification	<p>Recover viable mycobacterial specimens and monitor drug treatment.</p> <p>Identification and characterization of mycobacterial species.</p> <p>Isolation of MTBc for susceptibility testing.</p>	<p>More than 80% of MTBc and non-tubercular mycobacteria (NTM) are identified within 14 days of receipt.</p> <p>The first MTBc isolate from a new patient is communicated to submitter as <i>a critical value*</i>.</p> <p>Clinical specimens are allowed to incubate for 42 days (6 weeks) before being reported as culture negative.</p>
Antimycobacterial Susceptibility Testing (AST) for MTBc	<p>Determination of resistance or susceptibility of MTBc isolates to a variety of first line antimycobacterial drugs.</p>	<p>AST profiles are reported for more than 80% of isolates within 21 days of identification as MTBc. Resistances are communicated to submitter as <i>critical values*</i>.</p>
Submission to Collaborating Laboratories	<p>Specimens may be referred to collaborating labs for additional or confirmatory testing.</p> <p>In the case of rifampin resistance detected by NAAT, molecular detection of drug resistance (MDDR) analysis at the Centers for Disease Control and Prevention (CDC) can provide very rapid and broad initial susceptibility profiles.</p>	<p>Turnaround for results may vary depending on the collaborating lab.</p> <p>MDDR results are typically returned within seven days of package receipt at the CDC.</p>

Times to result will vary depending on quality of specimen collection, storage, transport, and media.


****Critical Values are time-sensitive results that must be communicated to the submitter by phone***




Mycobacteriology Lab Testing Sequence

Specimen Receipt at the Laboratory


(Turnaround times refer to specimens received Monday through before noon on Friday)




- Fluorescent smear (clinical specimen) or Kinyoun smear (referred culture) are reported as acid-fast bacilli seen or not; clinical specimens are processed and subcultured; incubation begins.
- New patient clinical specimens with acid-fast positive fluorescent smears are reported as *critical values*.



- NAA test performed, if requested, or if a new patient has an acid-fast positive smear from a clinical specimen. If rifampin resistance is detected, specimen will be forwarded to CDC for molecular detection of drug resistance (MDDR). Results are typically returned within seven days after package receipt at the CDC.
- All NAAT results, positive or negative, are reported as *critical values*.



- If mycobacteria are detected, species identification of MTBc or non-tuberculous mycobacteria (NTM) will be finished in about 80% of all specimens. Methods include both molecular techniques (PCR) and mass spectroscopy (MALDI-TOF).
- First isolates of MTBc for a new patient are reported as *critical values*.




- AST result for a new patient's first MTBc isolate completed in more than 80% of specimens.
- Any resistances are reported as *critical values*.

AT 42 DAYS – CLINICAL CULTURES COMPLETE INCUBATION. IF NO GROWTH IS SEEN IN ANY SUBCULTURE, SPECIMEN IS REPORTED AS “MYCOBACTERIA NOT FOUND”

If you still have additional questions regarding submitting mycobacterial specimens or testing, please contact us:

- Mycobacteriology Laboratory - phone (860) 920-6649 or fax (860) 920-6721.
- TB Epidemiology Program - phone (860) 509-7698 for courier information and other questions
- The Mycobacteriology Laboratory operates Monday through Friday, 7:30 am to 4:00 pm. The laboratory is closed weekends and on various federal and state holidays.



Appendix A

<p style="text-align: center;">◆ Submitter (REQUIRED)</p> <p>◆ LAB PROFILE Number:</p>	<p>CLINICAL TEST REQUISITION STATE OF CONNECTICUT Dr. Katherine A. Kelley State Public Health Laboratory 395 West Street, Rocky Hill, CT 06067 CLIA ID 07D0844555 Phone 860-920-6500 Form OL-9B Rev. 01/10/2024</p>	<p style="text-align: center;">ACCESSION LABEL FOR CTDPH LABORATORY USE ONLY</p>
<p>◆ DENOTES REQUIRED INFORMATION</p>		
<p>Section 1: Patient Information (Please Print Clearly)</p>		
<p>◆ Name (Last, First, M.I.) or Identifier:</p>		
<p>◆ Street Address:</p>		<p>◆ City, State, Zip:</p>
<p>◆ Date of Birth:</p>	<p>Gender: <input type="checkbox"/> Female <input type="checkbox"/> Male <input type="checkbox"/> Unknown</p>	<p>Home Phone:</p>
<p>Race (check all that apply): (◆ Race/Ethnicity Information is Required for Blood Lead)</p> <p><input type="checkbox"/> White <input type="checkbox"/> Black/African Amer. <input type="checkbox"/> Asian <input type="checkbox"/> Amer. Indian/Alaska Nat. <input type="checkbox"/> Nat. Hawaiian/Other Pacific Islander <input type="checkbox"/> Other <input type="checkbox"/> Unknown</p> <p>Ethnicity: <input type="checkbox"/> Hispanic <input type="checkbox"/> Non-Hispanic <input type="checkbox"/> Unknown</p>		
<p>◆ Ordering Healthcare Provider:</p>		<p>◆ Phone</p>
<p>Section 2: Specimen Information</p>		
<p>Specimen Storage (Prior to Delivery):</p> <p><input type="checkbox"/> Refrigerated (2-8°C) <input type="checkbox"/> Frozen (<-20°C) <input type="checkbox"/> Room Temperature</p>		<p>Specimen Receipt (CTDPH internal use only)</p> <div style="border: 1px solid black; width: 100px; height: 30px; display: flex; align-items: center; justify-content: center;"> °C </div> <p><input type="checkbox"/> Room Temperature <input type="checkbox"/> Refrigerated <input type="checkbox"/> Frozen</p>
<p>Specimen Transport/Delivery:</p> <p><input type="checkbox"/> Cold (Ice pack) <input type="checkbox"/> Frozen (Dry Ice) <input type="checkbox"/> Room Temperature</p>		
<p>Submitter Sample ID:</p>	<p>◆ Date Collected:</p>	<p>Time Collected: <input type="checkbox"/> AM <input type="checkbox"/> PM</p>
<p>◆ Specimen Source/Type:</p> <p><input type="checkbox"/> Axilla/groin <input type="checkbox"/> Blood (whole) <input type="checkbox"/> Bronchial Wash <input type="checkbox"/> Buccal cavity <input type="checkbox"/> Cervix <input type="checkbox"/> CSF <input type="checkbox"/> Lesion <input type="checkbox"/> Nasopharynx</p> <p><input type="checkbox"/> Oropharynx <input type="checkbox"/> Plasma <input type="checkbox"/> Rectum <input type="checkbox"/> Serum <input type="checkbox"/> Sputum <input type="checkbox"/> Stool <input type="checkbox"/> Urethra <input type="checkbox"/> Urine <input type="checkbox"/> Vagina</p> <p><input type="checkbox"/> Body Fluid, specify _____ <input type="checkbox"/> Tissue, specify _____</p> <p><input type="checkbox"/> Other, specify _____</p>		
<p>◆ Section 3: Select Testing Requested</p>		
<p>Bacteriology</p> <p><input type="checkbox"/> AFB Clinical Specimen (Mycobacteria Smear & Culture)</p> <p><input type="checkbox"/> AFB Referred Culture (Mycobacteria for Identification)</p> <p>Bacterial Isolate for Identification (Check one)</p> <p><input type="checkbox"/> Group A Streptococcus <input type="checkbox"/> Group B Streptococcus <input type="checkbox"/> <i>H. influenzae</i></p> <p><input type="checkbox"/> <i>L. monocytogenes</i> <input type="checkbox"/> Legionella <input type="checkbox"/> <i>N. meningitidis</i> <input type="checkbox"/> <i>S. pneumoniae</i></p> <p><input type="checkbox"/> Campylobacter <input type="checkbox"/> <i>E. coli</i> O157 <input type="checkbox"/> Salmonella <input type="checkbox"/> Shigella</p> <p><input type="checkbox"/> Shiga-toxin producing <i>E. coli</i> <input type="checkbox"/> Vibrio <input type="checkbox"/> Yersinia</p> <p><input type="checkbox"/> Other: _____</p> <p><input type="checkbox"/> Bioterrorism Agent Identification</p> <p>Specify agent: _____</p> <p><i>Bordetella pertussis</i> <input type="checkbox"/> Culture <input type="checkbox"/> DNA amplification</p> <p><input type="checkbox"/> Carbapenemase colonization screening (Rectal swab)</p> <p>Carbapenem resistant organism (Please attach susceptibility results)</p> <p><input type="checkbox"/> Fast Track (¹Epidemiology approval required)</p> <p><input type="checkbox"/> CRE (Enterobacterales, specify organism) _____</p> <p><input type="checkbox"/> CRAB (<i>Acinetobacter baumannii</i>) <input type="checkbox"/> CRPA (<i>Pseudomonas aeruginosa</i>)</p> <p><input type="checkbox"/> Enteric (Stool) Culture Suspect Organism: _____</p> <p><input type="checkbox"/> Shiga-toxin (+) Broth Culture</p>	<p>Serology/Virology/Sexually Transmitted Infections</p> <p>Arbovirus (Please select all that apply)</p> <p><input type="checkbox"/> Eastern Equine Encephalitis Virus IgM Antibody</p> <p><input type="checkbox"/> Powassan Virus IgM Antibody</p> <p><input type="checkbox"/> West Nile/St. Louis Virus IgM Antibody</p> <p><input type="checkbox"/> Chlamydia/ Gonorrhea Nucleic Acid Amplification Test</p> <p><input type="checkbox"/> Hepatitis A Virus PCR (¹Epidemiology approval required)</p> <p><input type="checkbox"/> Hepatitis B Surface Antibody</p> <p><input type="checkbox"/> Hepatitis B Surface Antigen</p> <p><input type="checkbox"/> Hepatitis C Testing</p> <p><input type="checkbox"/> Herpes Simplex IgG Antibody</p> <p><input type="checkbox"/> Herpes Simplex DNA amplification</p> <p><input type="checkbox"/> HIV-1/HIV-2 Antigen/Antibody</p> <p><input type="checkbox"/> HIV Viral Load</p> <p><input type="checkbox"/> Influenza/SARS-CoV-2 multiplex PCR</p> <p><input type="checkbox"/> Measles PCR</p> <p><input type="checkbox"/> MERS CoV PCR (Novel Coronavirus) (¹Epi Approval Required)</p> <p><input type="checkbox"/> Mumps PCR</p> <p><input type="checkbox"/> Non-Variola Orthopoxvirus PCR (R/O Monkeypox Virus)</p> <p><input type="checkbox"/> Norovirus PCR (¹Epidemiology approval required)</p> <p><input type="checkbox"/> QuantiFERON-TB Test (Specify ◆ Date AND Time Collected Above)</p> <p><input type="checkbox"/> Syphilis Screen (Serum)</p> <p><input type="checkbox"/> Syphilis VDRL (CSF)</p> <p><input type="checkbox"/> Respiratory Panel</p> <p><input type="checkbox"/> <i>Trichomonas vaginalis</i> NAAT (urine/vaginal Only)</p>	
<p>Blood Lead (Uninsured Patients ONLY) ◆ Race/Ethnicity Required</p> <p><input type="checkbox"/> Child Lead Screen (Capillary) <input type="checkbox"/> Confirmation (Venous)</p>		
<p>Mycology</p> <p><input type="checkbox"/> <i>Candida auris</i> identification (culture isolate)</p> <p><input type="checkbox"/> <i>Candida auris</i> screen</p> <p><input type="checkbox"/> Yeast ID/susceptibility testing (Blood <i>Candida</i> spp. Isolates ONLY)</p>		
<p>Parasitology</p> <p><input type="checkbox"/> Blood Parasite – Smear</p>		
<p>Comments</p> 		
<p>Test, Agent, or Disease Not Listed (Specify):</p>		
		<p>*Please provide:</p> <p>Symptoms _____</p> <p>Symptom onset date _____</p> <p>Travel history _____</p>
<p>¹DPH Epidemiology and Emerging Infections: (860)509-7994</p>		

Appendix B



Mycobacterium tuberculosis complex Nucleic Acid Amplification (NAA) Test Requisition

Katherine A. Kelley State Public Health Laboratory
395 West Street, Rocky Hill, CT 06067

For each clinical respiratory specimen where NAA testing is requested, complete this form, along with a Clinical Test Requisition, when submitting the specimen to the laboratory. Routine mycobacteria smear & culture will also be performed.

NAA testing will automatically be done on the first patient specimen submitted for routine mycobacteria smear & culture found to be **Acid-fast Bacilli (AFB) smear positive** by the CTDPH laboratory (the *M. tuberculosis* complex NAA Test Requisition is not required).

NAA Testing should **NOT** be ordered:

- When clinical suspicion is low (the positive predictive value of the test, the likelihood that the patient has tuberculosis when the test is positive, is low in such cases).
- To determine bacteriologic cure or to monitor response to antituberculous therapy

CTDPH TB Laboratory (Ph: 860-920-6649 / Fax: 860-920-6721)

CTDPH TB Control Program (860-509-7722)

<u>Submission Requirements</u>
<input type="checkbox"/> Clinical respiratory specimens (raw unprocessed): sputum, BAL, bronchial wash
<input type="checkbox"/> Patient has received no antituberculosis therapy, or less than three days of therapy at specimen collection.
<input type="checkbox"/> Specimens must be received by the laboratory within 10 days of collection.
<input type="checkbox"/> Test requests must be received within 7 calendar days of specimen receipt in the laboratory
<u>Submitter Information</u>
Authorized Submitter's Name: _____
Phone : _____ Fax: _____
<u>Patient Information</u>
Name: _____
Patient /Specimen ID #: _____ Date of Birth: _____
<u>Specimen Information</u>
Type / Source: <input type="checkbox"/> Sputum <input type="checkbox"/> Bronchoalveolar Lavage (BAL) <input type="checkbox"/> Bronchial Wash
Date Collected: _____ Other Information _____

Rev. 05/02/2024

Acknowledgements:


Rik Martinez

Christina Nishimura

Mark Harkins

Hongli Dong

Diane Charette



World TB Day



March 24