Pipette Verification	:			
Expected volume (ml)				
Replicate 1		Weight (g)	Weight (g)	Weight (g)
2 3 4 5				
6 7 8 9				
10	Replicate Average			
%Recovery Standard Deviation Coefficient of Variation	97.5-102.5% <5%			
Balance used:				
After reviewing the above	data it has been determine	ined that this pip	ette	
has met the acceptability criteria and is in proper working condition.				
has not met the acceptability criteria and is not in proper working condition.				
	Initial	l:	Date:	

Revision: 1

Effective: 8/20/2014

TX 14.2 Pipette Verification

Page 1 of 1

Note: Any pipettes which are marked as not in proper working condition will be immediately taken out of service and reported to the section supervisor.