| Chemical Preparedness and Response Revised 10-24-2023 | |
|---|--|
| Test description | Referral of blood and urine specimens to CDC for rapid toxic screen to assess the potential of a chemical exposure event; to provide surge capacity analysis, as needed, of human specimens for metabolites indicative of exposure to chemical agents. |
| Test Use | To aid in identifying a chemical terrorism event. |
| Test Department | Chemical Preparedness Phone: (860) 920-6716 or (860) 620-6717 |
| Methodology | ICP-MS, GC-MS |
| Availability | Daily, Monday-Friday |
| Specimen requirements | From each patient submit, Whole blood collected by venipuncture into 3 purple top tubes and 1 green or gray top tube Clean catch urine (40-60 ml) in a screw capped collection cup |
| | Blank tubes -2 empty and unopened purple top tubes, 2 empty and unopened green or gray top tubes, 2 empty and unopened urine cups |
| Collection | Collection kit containing tubes, submission form and collection instructions are |
| Kit/Container | located within the emergency departments in the state's sentinel hospitals |
| Collection | Standard venipuncture |
| Instructions | Clean catch urine collection |
| Specimen Handling & Transport | Transport to the laboratory immediately. |
| Unacceptable | Unlabeled specimens |
| conditions | Specimens that have leaked or containers that have broken in transit Improperly collected specimens |
| Requisition | Chemical Response Specimen Submission/ Chain of Custody Form |
| Form | |
| Required Information | Name and address of submitter. Two patient identifiers (Name, DOB, Acc.#, MRN), town of residence (city, state, zip), race/ethnicity, specimen type or source of collection, date collected, and test requested. Please ensure information on the requisition matches the specimen. |
| Additional Comments | The Chemical Preparedness and Response Laboratory is a partner in CDC's Chemical Laboratory Response Network (LRN-C) to develop surge capacity in response to a chemical exposure public health emergency. |