

*Approved by Director: Dr. Guy Vallaro***A. Purpose:**

To define symbols, abbreviations and acronyms that are commonly used in the latent print section.

B. Responsibility:

Latent Print Examiners

Abbreviation	Key
Δ	Change to
Ø	Identified
#	Number
?	Undetermined at this time
+	Positive
-	Negative
\	Right Sloped Loop
/	Left Sloped Loop
A/R	Administrative Review
ABP	Amino Acid based processing
AEL (laser)	Accessible Emission Level*
ALS	Alternate Light Source
AN, ANN	Annotated
AR	Administrative Review
ARCH ? P	Undetermined if Plain Arch
ARCH ? T	Undetermined if Tented Arch
AW	Accidental whorl
BBC	Black Backer Card
BC	Backer Card (BBC or WBC)
BI	Bi-chromatic
B/P or BP	Black Powder
BER	Blood Enhancement Reagent
BHL	Black Hinge Lifter
BLK	Black
BLK PDR	Black Powder
BLS	Blood Like Stain
BMAG	Black Magnetic Powder
BSA	By Submitting Agency
BP	Black Powder
CA	Superglue / Cyanoacrylate
CAE	Superglue / Cyanoacrylate Ester

LP SOP-05 Abbreviation & Symbol Key

Document ID: 1201

Revision: 11

Effective Date: 06/20/2023

Status: Published

Page 2 of 5

Approved by Director: Dr. Guy Vallaro

CAL	Calibrated
CHL	Clear Hinge Lifter
CM	Case Management
COMP	Comparison
CT-AFIS	Connecticut AFIS Database
CV	Comparison Value
CW (laser)	Continuous Wave*
CW	Central Pocket Loop Whorl
DEV	Developed
dH ₂ O	Distilled Water
DOB	Date of Birth
DPP	Dry powder processing
DSP	Dye stain processing
DW	Double Loop Whorl
ENF	Eliminations needed first
ENH	Enhanced
FA	Firearms
FB	Forensic Biology
F/B	Front & Back
FCR	Full Case Review
FF1	Forensic Filter range 532nm narrow band pass filter
FLR	Fluorescent
FP	Fingerprint
FRAGS	Fragments
FRD	Friction Ridge Detail
FRDD	Friction Ridge Detail - Documented
FT	Fingertip
GMAG	Grey Magnetic Powder
GP or G/P	Grey Powder
GRY	Grey
GYR	Green, Yellow, Red mark-up per SOP-LP-18
ID	Identification
IMP	Impression(s)
INC	Inconclusive
INSUFF	Insufficient
IQQ	Insufficient Quantity and Quality
IR	Infrared
IV	Identification Value
J (laser)	Joule*
JNT	Finger joint/phalanges
L	Left
L/L	Latent Lift
L/P or LP	Latent Print

LP SOP-05 Abbreviation & Symbol Key

Document ID: 1201

Revision: 11

Effective Date: 06/20/2023

Status: Published

Page 3 of 5

Approved by Director: Dr. Guy Vallaro

LAT REG	Latent Registered into Unsolved Latent Database
LFFS	Latent File Feature Search
LFIS	Latent File Image Search
LI	Left Index
LJ	Lower Joint
LL	Left Little
LM	Left Middle
LP	Latent Print
LPP	Left Palm Print
LR	Left Ring
LS/	Left Sloped
LSO (laser)	Laser Safety Officer*
LSP (laser)	Laser Safety Plan*
LT	Left Thumb
MCP	Major case prints
MPE (laser)	Maximum Permissible Exposure*
N.G.	Negative
N/A	Not Applicable
NAQ	Not AFIS Quality
NDVF	No Friction Ridge Detail of Value Found
NEATT	Not Examined At This Time
Neg	Negative
NG	Negative
NGI	Next Generation Identification
NHLR	No Hit Latent Registered
NHZ (laser)	Nominal Hazard Zone*
NI	Not Identifiable
NIN	Ninhydrin
NL	No Latents Found
NM	Nanometers
NO ANN	No annotation
NV	No Value
N/POR	Not porous
OCB	Orange Curved Barrier filter
OML	Oblique Magnified Lighting
ORIG	Original
ORO	Oil Red O
OV	Of Value
PC	Photocopy
PDR	Powder
PPI	Pixels per inch
POR	Porous
POS	Positive

LP SOP-05 Abbreviation & Symbol Key

Document ID: 1201

Revision: 11

Effective Date: 06/20/2023

Status: Published

Page 4 of 5

Approved by Director: Dr. Guy Vallaro

POSS	Possible
PP	Palm Print
PPC	Palm Print card
PW	Plain Whorl
QUAL	Quality
R or Rt	Right
R/B	Reddish/Brown
RC	Ridge Characteristics
RF or REF	Reference
Rev.	Reversed
RI	Right Index
RI	Rhode Island
RL	Right Little
RM	Right Middle
RPP	Right Palm Print
RPT	Report
RR	Right Ring
RS\	Right Sloped
RT	Right Thumb
RT	Right
SAF	Sample Accounted For (when using FB worksheets)
SBBHL	Sirchie Black Backer Hinge Lifter
SBP	Sebaceous based processing
sdH ₂ O	Sterile Distilled Water
SG	Superglue / Cyanoacrylate
SI	Superimposed
SID	State Identification Number
SPBI	State Police Bureau of Identification
SPR	Small Particle Reagent
SQQ	Sufficient Quantity and Quality
STHL	Sirchie Transparent Hinge Lifter
SUB	Submission
SWBHL	Sirchie White Backer Hinge Lifter
TIP	Finger tips
TL	Tape Lift
TL (laser)	Threshold Limit*
T/P	Test Print
T/R, T.R. or TR	Technical Review
THL	Transparent Hinge Lifter
TOT	Turned over To
TP	Test Print
TPC	Ten Print Card
TR	Technical Review

LP SOP-05 Abbreviation & Symbol Key

Document ID: 1201

Revision: 11

Effective Date: 06/20/2023

Status: Published

Page 5 of 5

Approved by Director: Dr. Guy Vallaro

TRANS	Transparent
UAL	Unknown Anatomical Location
UCN	Universal Control Number
ULD	Unsolved Latent Database
ULDB	Unsolved Latent Database
ULF	Unsolved Latent File
ULM	Unsolved Latent Match
UNK	Unknown
UV	Ultraviolet
VAL	Value
VE	Visual Examination
VER or VR	Verification or Verification Review
VIS	Visible
w/	With
W	Whorl
WBC	White Backer Card
WHL	White Hinge Lifter
WHT	White
WP	White Powder
WPP	Wet powder processing

*** denotes any abbreviations associated with the Laser and its use. Refer to LP SOP-28 and 29 for definition of abbreviations.**